



Neuza Salgado <sup>1,2</sup>, Mafalda Alexandra Silva <sup>1,3</sup>, Maria Eduardo Figueira <sup>4</sup>, Helena S. Costa <sup>1,3,\*</sup> and Tânia Gonçalves Albuquerque <sup>1,3</sup>

- <sup>1</sup> Research and Development Unit, Department of Food and Nutrition, National Institute of Health Dr. Ricardo Jorge, Avenida Padre Cruz, 1649-016 Lisbon, Portugal; tania.albuquerque@insa.min-saude.pt (T.G.A.)
- <sup>2</sup> Faculty of Pharmacy, University of Lisbon, Avenida Professor Gama Pinto, 1649-003 Lisbon, Portugal
   <sup>3</sup> REQUIMTE-LAQV/Faculty of Pharmacy, University of Porto, Rua Jorge Viterbo Ferreira, 228, 4050-313 Porto, Portugal
- <sup>4</sup> Research Institute for Medicines and Pharmaceutical Sciences (iMed.UL), Faculty of Pharmacy, University of Lisbon, Avenida Professor Gama Pinto, 1649-003 Lisbon, Portugal
- \* Correspondence: helena.costa@insa.min-saude.pt; Tel.: +351-217-519-200

Abstract: Oxalate is an antinutrient present in a wide range of foods, with plant products, especially green leafy vegetables, being the main sources of dietary oxalates. This compound has been largely associated with hyperoxaluria, kidney stone formation, and, in more severe cases, systematic oxalosis. Due to its impact on human health, it is extremely important to control the amount of oxalate present in foods, particularly for patients with kidney stone issues. In this review, a summary and discussion of the current knowledge on oxalate analysis, its extraction conditions, specific features of analytical methods, reported occurrence in foods, and its health implications are presented. In addition, a brief conclusion and further perspectives on whether high-oxalate foods are truly problematic and can be seen as health threats are shown.

Keywords: oxalate; foods; extraction conditions; analytical methods; occurrence; health implications

# 1. Introduction

Oxalate is a chemical compound that can form soluble and insoluble salts in water. This substance is present in a wide range of foods, with plant products being the main sources of dietary oxalates [1]. In plants, it plays a relevant role in various functions such as calcium homeostasis; pH regulation; plant growth, development and protection; photosynthesis; and detoxification of heavy metals [2,3].

According to the literature, various methods have been employed for the determination of oxalate in foods, including enzymatic assays [4–11], spectrofluorimetry [12], spectrophotometry, amperometry [9,13,14], electrochemical [15,16], capillary electrophoresis [17,18], titration [19–27], gas chromatography (GC) [28], and high-performance liquid chromatography (HPLC). HPLC is the most recently referenced method used for the determination of oxalates because of its high sensitivity, accuracy, versatility, and reliability, despite being expensive to purchase, repair, and maintain [29]. Conversely, despite its lower sensitivity, spectrophotometry is an inexpensive, rapid, and simpler method, requiring only one main instrument [30,31]. Accurate measurement of oxalate in foods is extremely dependent on its extraction, the first step in oxalate analysis [1]. Despite being two completely different methods, HPLC and spectrophotometry extraction conditions of total and soluble oxalates in foods are similar. Regarding analytical conditions, they have different and specific features since they are completely distinct procedures.

Concerning oxalate occurrence in foods, green leafy vegetables are particularly relevant, and some are considered high-oxalate foods. For example, published oxalate values are 329.6–2350 mg total oxalates/100 g fresh weight (FW) for spinach [32–39], 1235 mg total



Citation: Salgado, N.; Silva, M.A.; Figueira, M.E.; Costa, H.S.; Albuquerque, T.G. Oxalate in Foods: Extraction Conditions, Analytical Methods, Occurrence, and Health Implications. *Foods* **2023**, *12*, 3201. https://doi.org/10.3390/ foods12173201

Academic Editor: Mehdi Rahimmalek

Received: 31 July 2023 Revised: 19 August 2023 Accepted: 23 August 2023 Published: 25 August 2023



**Copyright:** © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). oxalates/100 g FW for rhubarb [1], 874 [40] and 1458.1 [1] mg total oxalates/100 g FW for swiss chard, 1079 mg soluble oxalates/100 g FW for sorrel [41], and 300.2–721.9 mg total oxalates/100 g FW for taro leaves [34].

Considering human health, oxalates have been a concern for a long time due to their antinutritive effects and potential nephrotoxicity [42,43]. As antinutrients, oxalates restrict the bioavailability of some nutrients since they can bind to minerals, reducing their absorption and use [3,44]. Potentially toxic soluble oxalates are delivered to the kidneys and can form calcium oxalate crystals there, which can lead to hyperoxaluria and kidney stone formation, also known as nephrolithiasis or urolithiasis [3,45–48]. In more serious cases, systemic oxalosis has been reported, a phenomenon in which calcium oxalate crystals deposit in various organs, tissues, and bones, when renal function declines and excess oxalate exists in the bloodstream [49,50].

This review has the purpose of gathering a considerable amount of information about oxalate, focusing on its extraction and analytical conditions and content in various foods measured by HPLC and spectrophotometry. Optimization of these parameters for oxalate determination in foods is very relevant to achieve reliable and accurate results considering the impact that this antinutrient could have on human health, especially on patients with kidney stone problems.

#### 2. Oxalates

Oxalate, or oxalic acid, is an antinutrient present, commonly in trace amounts, in fruits, nuts, cereals, fungi, vegetables, aromatic plants, and beverages, with plant-based products being the main sources of dietary oxalates [1]. However, some plants have high quantities of these compounds. In this matter, green leafy vegetables, such as spinach, Swiss chard, and rhubarb, are highlighted [3,32].

Oxalate can form soluble and insoluble salts in water. When binding with sodium, potassium, and ammonium ions, it forms soluble oxalates, whereas with calcium, iron, and magnesium it precipitates, forming insoluble compounds and making these minerals unavailable for absorption. Despite this fact, for example, zinc absorption and metabolism do not appear to be affected. In general, insoluble salts in water can be freely dissolved in acid [44,51,52]. Regarding health, one of the most important insoluble salts is calcium oxalate, having two hydration forms, monohydrates and dihydrates, which impact the shape of its crystals [1].

Depending on the pH of the cell sap, the liquid inside the vacuole of plants where oxalates are mostly found, they can present different chemical structures (Figure 1). On the one hand, when pH is 2, acid oxalate is the main oxalate. On the other hand, when pH is approximately 6, oxalate ion is the majority [51]. At the cytoplasmic pH of 7, oxalic acid also suffers deprotonation and exists as oxalate ion [2].



**Figure 1.** Chemical structure of (**a**) oxalic acid; (**b**) oxalate ion; (**c**) acid oxalate ion; and (**d**) oxalate salt, being  $M^{2+}$  a metallic cation.

Chemically, oxalic acid is characterized as a dicarboxylic organic acid with low molecular weight, high acidity ( $pK_{a1} = 1.25$ ,  $pK_{a2} = 4.27$ ), and chelating and reducing abilities. Therefore, in plants, it plays a relevant role in many biological processes such as calcium homeostasis; pH regulation; plant growth, development, and protection; photosynthesis; and detoxification of heavy metals [2,3]. However, when in excess in plants because of a metabolic disorder, this will promote impairment of its functions and, thus, reduction of

crop quality [2]. Many factors can influence oxalate accumulation in plants, such as growth, ripeness, variety, season, time of harvest, and cultivation conditions (e.g., use of nitrate fertilizer or soil conditions) [36,39,47,51].

The biosynthesis of oxalate in plants can result from different mechanisms, with glyoxylate, ascorbic acid, and oxaloacetate being the precursors of oxalate, an end product of their metabolisms. Therefore, there are three most-studied pathways: the glycolic acid/glyoxalic acid pathway, the ascorbic acid pathway, and the oxaloacetic acid pathway [2,51].

In addition to photosynthetic organisms, mammals can also produce oxalates in small amounts. In mammals, oxalate produced endogenously is a metabolite of ascorbate, hydroxyproline, glyoxylate, and glycine [3].

#### 3. Analysis of Oxalates in Foods

According to the literature, various methods have been employed for the determination of oxalate in foods, including enzymatic assays, spectrofluorimetry, spectrophotometry, amperometry, electrochemistry, capillary electrophoresis, titration, GC, and HPLC [53]. All of these methods have advantages and disadvantages like their high sensitivity and specificity but also high costs and time-consuming and complex sample handling [54]. For example, samples require an additional step (esterification) for GC analysis [28,55]. Even though there are many others, this manuscript will only focus on HPLC and UV–Vis spectrophotometry, specifically on their extraction and analytical conditions. HPLC is the most recently referenced method used for the determination of oxalates because of its high sensitivity, accuracy, versatility, and reliability, despite requiring equipment that is expensive to purchase, repair, and maintain [29]. Conversely, despite its lower sensitivity, spectrophotometry is an inexpensive, simple, rapid, and accurate technique, using only a spectrophotometer as the main instrument [30,31].

#### 3.1. Extraction Conditions

Accurate measurement of oxalate in foods is extremely dependent on its extraction, the first step in oxalate analysis. However, it seems to be difficult because of its extraction from plant tissue or its generation due to the oxidation of ascorbic acid during extraction [1].

Commonly, total oxalates (which include soluble and insoluble) are extracted with hydrochloric acid (HCl), whereas for soluble oxalates water is used [51]. However, a few papers describe different solutions, e.g., metaphosphoric acid [56], potassium phosphate buffer (pH 2.4) [57], or HCl with drops of octanol [58], for the extraction of total oxalates, and carbonate and sodium bicarbonate solution for the extraction of soluble oxalates [59,60]. Altunay et al. [30] also used HCl for the extraction of total oxalates and water for soluble oxalates; however, both extractions were conducted under ultrasonic power (300 W, 50 Hz).

The volume of the added extraction solution depends on the sample quantity, presenting a wide range of different values. Regarding other parameters, for HPLC analysis the hot extraction of total and soluble oxalates at 80 °C is more common, whereas for spectrophotometry 100 °C is more recurrent. For both methods, the most frequent time for extraction is 15 min. However, for HPLC, extraction times from 1–180 min (Table 1) are reported, whereas, for spectrophotometry, 15–1200 min is indicated (Table 2).

0	Matala	Sample		Extractio	on Condi	tions		Further Procedures	D (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	pН	Temp. (°C)	Time (min)	before Injection	Reference
	Rhubarb petioles	25	HCl (1 N)	200	-	100	15	Filtration and SPE (Sep-Pak cartridges) regenerated with 4 mL of methanol and 2 mL of water. The first 2 mL was discarded and the rest used for HPLC analysis.	[55]
	Carambola and spinach	100	HCl (3 N)	100	-	-	1	Filtration with Buchner funnel and two extractions by stirring each time with 50 mL of water or 3 N HCl. Concentration to 100 mL at 30 °C under reduced pressure (20 mmHg). SPE (C-18 Sep-Pak cartridge) pretreated with 2 mL of acetonitrile and 5 mL of water. The first 2 mL of eluate was discarded and the rest was filtered through a 0.45 $\mu$ m Millipore filter. Purification of some samples with ion exchange prior to HPLC analysis.	[35]
	Vegetables and spices	- 1			0.09			Centrifugation at 3000 rpm and filtration of 10 mL	[32,33]
	Vegetables, cereal grains, and legume seeds		HCl (2 N)	50	1	80	15	of the supernatant through a 0.45 mm cellulose acetate membrane.	[61]
T-1-1	Spinach	1	HCl (0.05 N)	-	-	-	-	-	[62]
Total			HCl (0.5 N)	5 (+5 deionized water)	-	100	20	Centrifugation at $12,000 \times g$ for 10 min. Supernatant removed and the volume increased to 20 mL with deionized water. Filtration with a filter with pore size 0.22 µm.	[39]
			HCl (25 %)			80			
			1101 (23 70)	_		100	30 and 180		
	Fruits, vegetables, beverages,					100		_	
	spices, herbs, nuts, cereals, algae, and mushrooms	2		4	-	21	15	Filtration of 1 mL of the solutions.	[1]
	algae, and mushrooms		HCl (2 N)				30	_	
						80	15	_	
							30		
	Cereal, cereal products, and plants of the Fabaceae, Convolvulaceae, and Malvaceae families	2	HCl (2 N)	4	-	21	15	Filtration.	[63,64]

**Table 1.** Extraction conditions for quantification of oxalates in foods by high-performance liquid chromatography.

0.14	<b>N</b> ( ( )	Sample		Extractio	n Condi	tions		Further Procedures	D (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	pН	Temp. (°C)	Time (min)	before Injection	References
	Spinach		HCl	40	0.93– 5.81	25–95	15–120	Cooled by standing in running cold water (11.5 °C) for 5 min. Volume made up to 100 mL with the extraction solution. Filtration of the sample solutions using a 0.45 $\mu$ m syringe filter into 1 mL glass vial.	[65]
	Taro leaves and corms and Indian vegetables	3 or 0.5	HCl (2 N)	50	0.7	80	15	Centrifugation at 3000 rpm and filtration of 10 mL of the supernatant through a 0.45 mm cellulose acetate membrane.	[34]
	Spinach	0.3	HCl (0.2 N)	10		80	15	Centrifugation at $2889 \times g$ for 15 min. Filtration of the supernatant through a 0.45 mm cellulose nitrate filter.	[66]
	Green juices with spinach and other vegetables and fruits	_	HCI (0.2 IN)	40	-	80	20	The extracts were allowed to cool and then made up to 100 mL, in a volumetric flask, with HCl (0.2 N).	[67]
	Korean vegetables	- 5	HCl (2 N)	40 (+20 after homogenization)	-	80	15	Filtration with filter paper. Dilution with deionized water (90 mL). Filtration with 0.45 µm regenerated cellulose microfilters.	[40]
	Pakistani vegetables and beans	5 (mixed with 5 mL of distilled water)	HCl (2 N)	50	-	-	-	Centrifugation at 5000 rpm for 20 min and the supernatants were transferred to 100 mL volumetric flasks and made up to final volume with distilled water.	[36]
	Taro corms	0.5	HCl (2 N)	20	-	80	15	The extracts were cooled at room temperature and made up to 50 mL with HCl (2 N). Centrifugation at 2889 rpm for 15 min. Separation and filtration of the supernatant through a 0.45 mm cellulose nitrate membrane.	[68]
	Mexican vegetables	_	HPO <sub>3</sub> (45 g/L)	25	-	-	15	Centrifugation (30 min, $8960 \times g$ ) and the supernatant made up to 25 mL with metaphosphoric acid.	[56]
	Ethiopian collard greens and mango	4	HCl (2 N)	50	-	-	50 (at 250 rpm)	The solution was removed from the shaker and 50 mL HPLC grade $H_2O$ was added. Filtration through 0.45 $\mu$ m syringe filter. The filtrate was transferred into a 2 mL vial.	[58]
	Various foods	2 to 3	HCl (0.2 N)	-	-	60	60	Centrifugation (15,000 $\times$ g; 5 min) and filtration through 25 mm diameter 0.2 $\mu$ m PTFE filter.	[17]

0.14	N ( )	Sample		Extraction	n Conditi	ions		Further Procedures	<b>D</b> (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	pН	Temp. (°C)	Time (min)	before Injection	References
	Spinach and kale	-	Potassium phosphate buffer (0.2 M)	-	2.4	-	-	-	[57]
	Carambola, spinach, spinach products, and New Zealand spinach	100	Water	100	-	-	1	Filtration with Buchner funnel and two extractions by stirring each time with 50 mL of water or 3 N HCl. Concentration to 100 mL at 30 °C under reduced pressure (20 mmHg). SPE (C-18 Sep-Pak cartridge) pretreated with 2 mL of acetonitrile and 5 mL of water. The first 2 mL of eluate was discarded and the rest was filtered through a 0.45 $\mu$ m Millipore filter. Purification of some samples with ion exchange prior to HPLC analysis.	[35,69]
	Spinach		Carbonate (1.8 mM) and sodium	50	-	-	-	-	[60]
	Fresh vegetables		bicarbonate (1.7 mM) solution		-	25	30	-	[59]
Soluble	Taro corms	0.5	Deionized water	20	-	80	15	The extracts were cooled at room temperature and made up to 50 mL with deionized water. Centrifugation at 2889 rpm for 15 min. Separation and filtration of the supernatant through a 0.45 mm cellulose nitrate membrane.	[68]
	Mexican vegetables		Distilled water	25	-	-	15	Centrifugation (30 min, $8960 \times g$ ) and the supernatant was made up to 25 mL with metaphosphoric acid.	[56]
	Vegetables and spices				6.90			Centrifugation (3000 rpm) and filtration of 10 mL of	[32,33]
	Vegetables, cereal grains, and legume seeds	1	Distilled water	50	5.50	80	15	the supernatant through a 0.45 mm cellulose acetate membrane.	[61]
	Samples of herbal, green, oolong, and black teas	-	Distilled water	50	6.90	80	15	Each sample of tea was filtered through a 0.45 mm cellulose acetate membrane syringe filter.	[70]
	Taro leaves and corms and Indian vegetables	3 or 0.5	Distilled water	50	6.50	80	15	Centrifugation at 3000 rpm and filtration of 10 mL of the supernatant through a 0.45 mm cellulose acetate membrane.	[34]

0.17		Sample		Extractio	n Condi	tions		Further Procedures	<b>D</b> (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	pН	Temp. (°C)	Time (min)	before Injection	Reference
						21	15		
	Fruits, vegetables, beverages,					21	30	_	
	spices, herbs, nuts, cereals,	2	Distilled water	4	-	20	15	<ul> <li>Filtration and acidification by adding HCl (50</li> <li>μL/mL, 2 N) to stabilize ascorbic acid.</li> </ul>	[1]
	algae, and mushrooms		water			80	30		
					-	21	15	_	
	Green and black tea		Water	-	-	100	6	Filtration of the obtained extract through a disposable 0.45 $\mu$ m filter. Dilution 10 times.	[71]
	Cereal, cereal products, and plants of the Fabaceae, Convolvulaceae, and Malvaceae families	2	Distilled water	4	-	21	15	Filtration and acidification by adding HCl (50 $\mu$ L/mL, 2 N) to stabilize ascorbic acid.	[63,64]
	Korean vegetables		Water	40 (+20 after homogenization)	-	80	15	Filtration with filter paper. Dilution with deionized water (90 mL). Filtration with 0.45 µm regenerated cellulose microfilters.	[40]
	Green juices with spinach and other vegetables and fruits	5	Nanopure II water	40	-	80	20	The extracts were allowed to cool and then made up to 100 mL, in a volumetric flask, with Nanopure II water.	[67]
	Pakistani vegetables and beans	5 (mixed with 5 mL of distilled water)	Distilled water	50	-	-	-	Centrifugation at 5000 rpm for 20 min and the supernatants were transferred to 100 mL volumetric flasks and made up to final volume with distilled water.	[36]
	Spinach	0.3	Nanopure water	40	-	80	15	The extracts were allowed to cool and then transferred into 100 mL volumetric flasks and made up to final volume. Centrifugation at $2889 \times g$ for 15 min. The supernatant was filtered through a 0.45 mm cellulose nitrate filter.	[66]
	Alcoholic and non-alcoholic beverages	1.75 to 7	Water	200 or 150	-	70	5	Filtration of 4 mL of samples and acidification with 50 μL HCl (2 N).	[72]

0.17	Matrix	Sample		Extractio	n Condi	tions	Further Procedures	<b>D</b> (	
Oxalates		Amount (g)	Solution	Volume (mL)	pН	Temp. (°C)	Time (min)	before Injection	References
	Spinach	-	Deionized water	10 times the volume of fresh weight	_	RT	5	Filtration through filtrate paper (ADVANTEC no. 3; Advantec, Tokyo) on ice. Further filtration using the spin column Ultra-free-MC 0.45 $\mu$ m PTFE membrane (Millipore, Billerica, Mass.) and centrifugation at 5000 × g for 60 min at 4 °C. Storage of samples at -20 °C. Frozen filtrates were thawed on ice, diluted to 10× with deionized water, and transferred to a disposable vial (S/T micro vial; Tomsic, Tokyo).	[73]
	Spinach	1	Deionized water	5 (+5 after cooling)	-	100	20	Centrifugation at $12,000 \times g$ for 10 min. The supernatant was removed and the volume made up to 20 mL with deionized water. Filtration with a filter with pore size 0.22 µm.	[39]

RT—room temperature; Temp.—temperature; "-"—information not available.

 Table 2. Extraction conditions for quantification of oxalates in foods by spectrophotometry.

0.11	N ( )	Sample		Extraction Condit	ions		- Further Procedures before Analysis	References
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	Temp. (°C)	Time (min)	- Further Frocedures before Analysis	Kererences
	Cress, dill, parsley, cauliflower, broccoli, celery, black cabbage, red radish, lettuce, and leek	3	HCl (0.2 N)	30	60	15	The mixture of the sample with HCl was degassed and digested under ultrasonic power (300 W, 50 Hz). It was allowed to cool and then filtered by using a membrane filter of $0.45 \ \mu m$ into a 100 mL volumetric flask. The final volume was diluted to 100 mL with ultrapure water before analysis.	[30]
Total	Ethiopian collard greens, cabbage, lettuce, beetroot, pineapple, and mango	1	(a) Water, (b) HCl (6 N), (c) octanol	(a) 150, (b) 27.5, (c) 2 drops	100	25	The mixture was cooled, transferred to a 250 mL volumetric flask, and the volume completed. Filtration through Whatman 541 filter paper. Evaporation of 10 mL of this filtrate at 40–45 °C in a vacuum oven and redissolution in 10 mL of 0.01 M H <sub>2</sub> SO <sub>4</sub> .	[58]
	Spinach	0.50	HCl (2 N)	20	-	-	Centrifugation at 2500 rpm for 6 min and the supernatant was filtered through Whatman No. 1 paper. The residue retained by the filter was treated twice with 10 mL of HCl (2 N) and the combined filtrates were	[37]
		5	HCl (2 N)	20	-	-	diluted to 250 mL with water.	[38]

0.1/		Sample		Extraction Condition	ons			D (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	Temp. (°C)	Time (min)	- Further Procedures before Analysis	References
	Spinach and mushroom	5	HCl (2 N)	20	-	-	Centrifugation of the suspension at 2500 rpm for 6 min. Filtration of the supernatant through filter paper (Whatman No. 1). The residue retained by the filter was treated twice with 10 mL of HCl (2 N) and then filtered. The combined filtrates were mixed and diluted to 100 mL with water.	[74]
	Spinach and mushroom	2.5 or 25	Water	-	100	20	The suspension, previously cooled, was centrifuged and then filtered through filter paper (Whatman No. 1). Dilution of the filtrate to 1000 cm <sup>3</sup> and 2.0 cm <sup>3</sup> of the sample solution was used in the proposed method.	[54]
	Cress, dill, parsley, cauliflower, broccoli, celery, black cabbage, red radish, lettuce, and leek	3	Ultrapure water	30	60	15	The mixture of the sample with water was degassed and digested under ultrasonic power (300 W, 50 Hz). It was allowed to cool and then filtered by using a membrane filter of 0.45 µm into a 100 mL volumetric flask. The final volume was diluted to 100 mL with ultrapure water before analysis.	[30]
	Sorrel, spinach, parsley, ginger, and black pepper	10	Distilled water	100	100	15	Filtration after cooling. The solution pH was adjusted to 3.0 and the final solution was diluted to 100 mL. Dilution five times.	[41]
Soluble	Tap water	50	-	-	100	20	-	[75]
	Spinach and mushroom	3 or 15	Water	Necessary to dilute to 100 mL in a calibrated flask	100	45	After cooling, the filtration of the suspension was carried out through Whatman No. 1 filter paper. Dilution to 250 mL. Adjustment of the pH to about 10 by dropwise addition of NaOH (0.1 N). The solution was centrifuged (MS-3400 centrifuge; Cole-Parmer, St. Neots, UK) at $1492 \times g$ for 5 min. After the neutralization of the solution with HCl (0.1 N), it was diluted in a 100 mL volumetric flask.	[76]
	Vegetables, beverages, and fruits	20	Distilled water	60 (+HCl to adjust to pH 2–3)	50	15	After cooling to room temperature and adjusting to pH 3.0 with potassium hydroxide, the mixtures were transferred to 100 mL volumetric flasks and they were made up with distilled water. Filtration through membrane filters (0.45 µm).	[77]

0.14	N ( )	Sample		Extraction Condition	ons		Furthern Duran Junior Instance Arrists	D (
Oxalates	Matrix	Amount (g)	Solution	Volume (mL)	Temp. (°C)	Time (min)	- Further Procedures before Analysis	References
	Spinach and mushrooms	3 or 15	Water	-	100	45	The suspension was filtered twice through filter paper (Whatman No. 1) and the filtrate was diluted to 250 mL. Adjustment of the pH to about 10 by dropwise addition of 0.10 N sodium hydroxide solution. Centrifugation at 2000 rpm for 5 min. The resulting solution was decanted, neutralized with HCl (0.10 N), and diluted in a 100 mL volumetric flask. Then, 1.0 mL of the solution obtained was used for the proposed method.	[78]
-	Spinach, mushroom, and kidney bean	5 to 10	Water	250	-	10	Centrifugation at 2500 rpm for 5 min. Filtration of the supernatant dryly through Whatman No. 1 paper.	[79]
	Spinach	5	Water	-	100	30	Filtration through a filter paper, after being cooled. The filtrate was diluted to 250 mL and 5 mL of this solution was used to determine oxalate.	[80]
	Spinach	5	Water	Necessary to dilute to 50 mL in a volumetric flask	-	1200	Centrifugation of the suspension at 2500 rpm for 10 min. Filtration of the supernatant through Whatman No. 1 paper.	[81]
·	Spinach	5	Water	Necessary to dilute to 50 mL in a volumetric flask	-	-	Centrifugation of the suspension at 2500 rpm for 5 min. Filtration of the supernatant through Whatman No. 1 paper.	[31]
-	Beetroot, spinach, and mushroom	50	Water	-	100	50	The mixture was cooled and filtered through a membrane filter.	[82]

"-"—information not available.

Additionally, some studies have shown the influence of modifying these extraction conditions. Hönow et al. [1] concluded that increasing extraction time from 30 to 180 min for the extraction of total oxalates resulted in oxalate generation and the increase was higher when treated at 100 °C (reflux) than at 80 °C (water bath). For soluble oxalates, results increased significantly after extraction at 80 °C compared to extraction at 21 °C and there were no significant differences between extractions of 15 or 30 min, proposing that soluble oxalate should be extracted with distilled water for 15 min at room temperature. In contrast, other authors suggested that room temperature might not be enough for the complete extraction of oxalates, leading to an underestimation [5]. This theory is also supported by Kusuma et al. [65] who considered that extraction temperatures above 65 °C are required for efficient extraction of total oxalates and pH should be at least 1. Therefore, the ideal temperature of oxalate extraction remains a controversial question because it can lead to oxalate generation due to in vitro conversion from precursors or failure to dissolve all pre-existing insoluble oxalates [1,5]. Concerning time, 15 min was considered to be the minimum for extraction and increasing it has not been associated with any advantage [65].

Frequently, further procedures after extraction mainly consist of filtration and centrifugation. However, some authors utilize more specific techniques. For example, prior to HPLC analysis, there have been reports of the use of solid-phase extraction (SPE) [35,55,69], purification with ion exchange [35], concentration [35,69], and, specifically for soluble oxalates, acidification with HCl to stabilize ascorbic acid which can be present and cause oxaloneogenesis at pH above 5, resulting in an overestimation [1,63,64,72]. For spectrophotometric methods, evaporation and subsequent redissolution for total oxalate analysis have been mentioned [58] and, once again, adjustment of pH, which is relevant for soluble oxalates [41,76,78].

Insoluble oxalates are also extracted during treatment with HCl, so their content is always calculated by the difference between total and soluble oxalates as Holloway et al. [83] suggested.

### 3.2. Methods' Conditions

Despite having similar extraction conditions, chromatographic and spectrophotometric methods have different and specific features (Tables 3 and 4).

Matrix	Analytes	Chromatographic Conditions	Validation Data	References
Rhubarb petioles	Total oxalates	Column: LiChrosorb RP-8 (250 $\times$ 4.6 mm; 10 $\mu$ m particle size) Guard column: LC C-18 guard column Detector ( $\lambda$ , nm): UV (220) Mobile phase: 0.5% KH <sub>2</sub> PO <sub>4</sub> and 0.005 M TBA buffered at pH 2.00 with orthophosphoric acid Type of elution: Isocratic Injection volume ( $\mu$ L): 5 Flow rate (mL/min): 2 Column temp: (°C): - Run time (min): 7	Linearity range: 0–0.40 mg/mL Determination coefficient (r <sup>2</sup> ): - LOD: < 0.001 mg/mL LOQ: -	[55]
Carambola and spinach	Total oxalates	Column: Zorbax amine $(250 \times 4.6 \text{ mm}; 7 \mu \text{m} \text{ particle size})$ Guard column: Amino guard column (10 $\mu \text{m} \text{ particle size})$ Detector ( $\lambda$ , nm): UV (206) Mobile phase: Buffer solution of aqueous NaH <sub>2</sub> PO <sub>4</sub> (0.15 M) at pH 2.4 Type of elution: Isocratic Injection volume ( $\mu$ L): 20 Flow rate (mL/min): 1.1 Column temp: (°C): - Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[35]
-	Organic acids (e.g., oxalic acid)	Columns: (1) Radial compression column with C-18 or C-8 functionality (100 mm; 10 $\mu$ m particle size); (2) C-8 analytical columns of two manufacturers (250 × 4.6 mm, 6 and 10 $\mu$ m particle size); (3) propylamine anion exchange column (25 cm × 4.6 mm; 6 $\mu$ m particle size) (4) diethylaminoethyl (DEAE) anion exchange column (250 × 4.6 mm; 5 $\mu$ m particle size); and (5) Hamilton HA-X8.00 column with strong anion exchange resin in the sulfate form (255 × 5 mm; 7–10- $\mu$ m particle size) Guard column: - Detector ( $\lambda$ , nm): RI and UV (254 and 206) Mobile phase: (1–2) 2% NH4H2PO4 adjusted to pH 2.4 with phosphoric acid; (3) 0.15 M NaH2PO4 (pH = 4.2); (4) 0.30 M NH4H2PO4 adjusted to pH 6.5 with concentrated ammonium hydroxide; (5) 0.5 M (NH4)2SO4 adjusted to pH 7.25 with ammonium hydroxide; 0.3 M (NH4)2SO4 containing 10% methanol and 0.1–1.5 M MgSO4 solutions (for the gradient analysis) Type of elution: (1–5) Isocratic and (5) linear gradient Injection volume ( $\mu$ L): 20 Flow rate (mL/min): 1; 1.5 and 2 Column temp: (°C): (columns 1 to 4) RT and (column 5) 80 Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[84]

# **Table 3.** High-performance liquid chromatography conditions for the analysis of oxalates in foods.

Matrix Analytes **Chromatographic Conditions** Validation Data References Column: IonPac AS4A (250  $\times$  4 mm; 15  $\mu$ m particle size) Guard column: IonPac AG4A Detector ( $\lambda$ , nm): Conductivity Mobile phase:  $Na_2CO_3$  (1.8 mM) and Linearity range: -NaHCO<sub>3</sub> (1.7 mM) solution Determination coefficient (r<sup>2</sup>): -[60] Spinach Soluble oxalates Type of elution: Isocratic LOD: -Injection volume (µL): -LOO: -Flow rate (mL/min): -Column temp: (°C): -Run time (min): -Column: IonPac AS4A (250  $\times$  4 mm; 15  $\mu$ m particle size) Guard column: IonPac AG4A Detector ( $\lambda$ , nm): Conductivity Mobile phase: Na<sub>2</sub>CO<sub>3</sub> (1.8 mM) and Linearity range: -NaHCO<sub>3</sub> (1.7 mM) solution Determination coefficient (r<sup>2</sup>): -[59] Vegetables Soluble oxalates Type of elution: Isocratic LOD: -Injection volume (µL): -LOQ: -Flow rate (mL/min): 2 Column temp: (°C): -Run time (min): -Column: Dionex IonPac AS4A-SC ( $250 \times 4$  mm; 13 µm particle size) Guard column: -Detector  $(\lambda, nm)$ : -Linearity range: -Mobile phase: -Determination coefficient (r<sup>2</sup>): -Spinach Total oxalates Type of elution: -[62] LOD: -Injection volume (µL): -LOQ: -Flow rate (mL/min): -Column temp: (°C): -Run time (min): -Column: Bio-Rad Aminex ion exclusion HPX-87H ( $300 \times 7.8$  mm; 9 µm particle size) Guard column: Aminex Cation-H guard column Detector ( $\lambda$ , nm): UV–Vis (210) Linearity range: 0.01–0.2 mg/mL Mobile phase:  $H_2SO_4$  (0.0125 M) filtered through a 0.45 µm membrane and degassed Spinach, swiss chard, Determination coefficient  $(r^2)$ : 0.999 using vacuum (total) and 0.986 (soluble) broccoli, carrot, parsnip, Total and soluble oxalates [32] Type of elution: Isocratic LOD: rhubarb stalks, and beetroot Injection volume ( $\mu$ L): 5 LOQ: -Flow rate (mL/min): 0.5 (0.1 prior to use and in between sample sets) Column temp: (°C): RT Run time (min): -

Matrix	Analytes	Chromatographic Conditions	Validation Data	References
Various foods	Soluble oxalates	$ \begin{array}{c} \mbox{Column: Alltech All-Sep anion exchange column (100 \times 4.6 mm; 7  \mu m \mbox{ particle size}) \\ \mbox{Guard column: -} \\ \mbox{Detector } (\lambda, nm): \mbox{Conductivity} \\ \mbox{Mobile phase: } Na_2 CO_3 \ (0.9 \mbox{ mM}) \mbox{ and } \\ \mbox{NaHCO}_3 \ (0.85 \mbox{ mM}) \\ \mbox{Type of elution: Isocratic} \\ \mbox{Injection volume } (\mu L): - \\ \mbox{Flow rate } (mL/min): 1.2 \\ \mbox{Column temp: } (^{\circ}C): - \\ \mbox{Run time (min): -} \end{array} $	Linearity range: $0.5 \times 10^{-9}$ – $1 \times 10^{-8}$ mg/mL Determination coefficient (r <sup>2</sup> ): 0.990 LOD: - LOQ: 0.2 mg/100 g	[17]
Samples of herbal, green, oolong, and black teas	Soluble oxalates	Column: Bio-Rad Aminex ion exclusion HPX-87H ( $300 \times 7.8 \text{ mm}; 9 \mu \text{m}$ particle size) Guard column: - Detector ( $\lambda$ , nm): UV (210) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (0.0125 M) Type of elution: Isocratic Injection volume ( $\mu$ L): 5 Flow rate (mL/min): 0.5 Column temp: (°C): 50 Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[70]
Fruits, vegetables, beverages, spices, herbs, nuts, cereals, algae, and mushrooms	Total and soluble oxalates	$ \begin{array}{c} \mbox{Column: Dionex IonPac AS4A anion exchange column (250 \times 4 mm; 15  \mu m \\  particle size) \\  Guard column: - \\  Detector (\lambda, nm): Amperometric \\  Mobile phase: 2.0 g EDTA/L distilled water adjusted to pH 5.0 by adding 15 mL of \\  0.3\%  NaOH  suprapur \\  Type  of  elution:  Isocratic \\  Injection  volume  (\mu L): - \\  Flow  rate  (mL/min): - \\  Column  temp: (^{\circ}C): - \\  Run  time  (min): - \\ \end{array} $	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: 0.68 µM LOQ: -	[1]
Spinach	Soluble oxalates	Column: Mightysil RP-C18 Aqua column (250 × 4.6 mm; 5 μm particle size) Guard column: - Detector (λ, nm): UV (210) Mobile phase: Tetrabutylammonium chloride (5 mM) in phosphate ammonium buffer (pH 6.8) Type of elution: Isocratic Injection volume (μL): - Flow rate (mL/min): 1 Column temp: (°C): 30 Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[73]

Matrix	Analytes	Chromatographic Conditions	Validation Data	References
Cereal and cereal products	Total and soluble oxalates	Column: Dionex IonPac AS4A anion exchange column ( $250 \times 4 \text{ mm}$ ; 15 µm particle size) Guard column: - Detector ( $\lambda$ , nm): Amperometric Mobile phase: 2.0 g of EDTA/L of distilled water adjusted to pH 5.0 by adding 15 µL of 0.3 N NaOH Type of elution: Isocratic Injection volume (µL): - Flow rate (mL/min): - Column temp: (°C): - Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[64]
Thai vegetables, cereal grains, and legume seeds	Total and soluble oxalates	Column: Bio-Rad Aminex ion exclusion HPX-87H ( $300 \times 7.8 \text{ mm}$ ; 9 µm particle size) Guard column: - Detector ( $\lambda$ , nm): UV (210) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (0.0125 M) Type of elution: Isocratic Injection volume (µL): - Flow rate (mL/min): - Column temp: (°C): - Run time (min): -	Linearity range: 0.1–0.5 mg/mL Determination coefficient (r <sup>2</sup> ): 0.997 LOD: - LOQ: 3 mg/100 g	[61]
Korean vegetables	Total and soluble oxalates		Linearity range: 0.0168–1.3131 mg/mL Determination coefficient (r <sup>2</sup> ): 0.9995 LOD: - LOQ: -	[40]
Taro leaves and corms and Indian vegetables	Total and soluble oxalates	$ \begin{array}{c} \mbox{Column: Bio-Rad Aminex ion exclusion HPX-87H (300 \times 7.8 mm; 9 \mbox{ \mum particle size}) \\ \mbox{Guard column: Aminex Cation-H guard column} \\ \mbox{Detector } (\lambda, nm): UV-Vis (210) \\ \mbox{Mobile phase: H}_2SO_4 (0.0125 \mbox{ M}) filtered through a 0.45 \mbox{ \mum membrane and degassed} \\ \mbox{using vacuum} \\ \mbox{Type of elution: Isocratic} \\ \mbox{Injection volume } (\mu L): 5 \\ \mbox{Flow rate (mL/min): 0.5 (0.1 prior to use and in between sample sets)} \\ \mbox{Column temp: (°C): RT} \\ \mbox{Run time (min): -} \end{array} $	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: 5 mg/100 g DM LOQ: -	[34]

Matrix	Analytes	<b>Chromatographic Conditions</b>	Validation Data	References
Pakistani vegetables and beans	Total and soluble oxalates	Column: Supelco reversed-phase column ( $250 \times 4.6 \text{ mm}$ ; 5 µm particle size) Guard column: - Detector ( $\lambda$ , nm): UV (210) Mobile phase: 0.25% dihydrogenate phosphate and 0.0025 M tetrabutylammonium hydrogen sulfate buffered at pH 2.0 with ortho-phosphoric acid Type of elution: Isocratic Injection volume (µL): 5 Flow rate (mL/min): 1 Column temp: (°C): - Run time (min): -	Linearity range: 1 × 10 <sup>-3</sup> –0.04 mg/mL Determination coefficient (r <sup>2</sup> ): 0.9773 LOD: - LOQ: -	[36]
Spinach	Total and soluble oxalates	Column: Rezex ROA ion exclusion organic acid column $(300 \times 7.8 \text{ mm}; 8 \mu \text{m} \text{particle size})$ Guard column: Bio-Rad cation-H guard column Detector ( $\lambda$ , nm): UV-Vis (210) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (25 mM) Type of elution: Isocratic Injection volume ( $\mu$ L): 20 Flow rate (mL/min): 0.6 Column temp: (°C): 25 Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[66]
Green and black tea	Soluble oxalates	Column: Shodex IC SI-90 anion exchange column $(250 \times 4 \text{ mm}; 9 \mu\text{m} \text{ particle size})$ filled with KanK-ASt $(120 \times 5 \text{ mm}, 14 \mu\text{m} \text{ particle size})$ Guard column: - Detector ( $\lambda$ , nm): Conductivity Mobile phase: Na <sub>2</sub> CO <sub>3</sub> (1.9 mM) and NaHCO <sub>3</sub> (2.4 mM) solution Type of elution: Isocratic Injection volume ( $\mu$ L): 20 Flow rate (mL/min): 1.5 Column temp: (°C): 33 Run time (min): -	Linearity range: $1 \times 10^{-4}$ -0.02 mg/mL Determination coefficient (r <sup>2</sup> ): 0.9998 LOD: $3 \times 10^{-5}$ mg/mL LOQ: $1 \times 10^{-4}$ mg/mL	[71]
Spinach	Total and soluble oxalates	Column: Hypersil C18 column ( $250 \times 4.6 \text{ mm}$ ; 5 µm particle size) Guard column: - Detector ( $\lambda$ , nm): UV (220) Mobile phase: Aqueous solution containing 0.5% KH <sub>2</sub> PO <sub>4</sub> and 0.5 mM tetra- <i>n</i> -butyl ammonium hydrogen sulfate (pH 2.0) degassed with an ultrasonic generator for 20 min Type of elution: Isocratic Injection volume (µL): 5 Flow rate (mL/min): 0.5 Column temp: (°C): 40 Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[39]

Matrix	Analytes	Chromatographic Conditions	Validation Data	References
Green juices with spinach and other vegetables and fruits	Total and soluble oxalates	Column: Rezex ROA ion exclusion organic acid column (300 × 7.8 mm; 8 μm particle size) Guard column: Bio-Rad cation-H guard column Detector (λ, nm): UV- <i>Vis</i> (210) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (25 mM) Type of elution: Isocratic Injection volume (μL): 20 Flow rate (mL/min): 0.6 Column temp: (°C): - Run time (min): -	Linearity range: 0.01–0.25 mg/mL Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[67]
Spinach	Total oxalates	Column: Rezex ROA ion exclusion organic acid column (300 × 7.8 mm; 8 μm particle size) Guard column: Bio-Rad cation-H guard column Detector (λ, nm): UV- <i>Vis</i> (210) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (25 mM) Type of elution: Isocratic Injection volume (μL): 20 Flow rate (mL/min): 0.6 Column temp: (°C): 25 Run time (min): -	Linearity range: 0.01–0.25 mg/mL Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[65]
Alcoholic and non-alcoholic beverages	Soluble oxalates	Column: Dionex IonPac AS4A anion exchange column (250 × 4 mm; 15 μm particle size) Guard column: - Detector (λ, nm): Amperometric Mobile phase: Aqueous EDTA solution (2.0 g/L) adjusted to pH 5.0 with 0.3 M NaOH Type of elution: Isocratic Injection volume (μL): - Flow rate (mL/min): 0.6 Column temp: (°C): - Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[72]

Matrix	Analytes	Chromatographic Conditions	Validation Data	References
Vegetables and fruits Total oxalates		Column: Agilent Poroshell-C18 (250 $\times$ 4.6 mm; 2.7 $\mu$ m particle size) Guard column: - Detector ( $\lambda$ , nm): UV (210) Mobile phase: 50 mM KH <sub>2</sub> PO <sub>4</sub> , H <sub>3</sub> PO <sub>4</sub> (pH 2.8) Type of elution: Isocratic Injection volume ( $\mu$ L): 5 Flow rate (mL/min): 1 Column temp: (°C): 50 Run time (min): 25	Linearity range: 0.1–1.0 mg/mL Determination coefficient (r <sup>2</sup> ): -	(=0)
	Iotal oxalates	Column: Agilent Poroshell-C18 (100 $\times$ 2.1 mm; 2.7 µm particle size) Guard column: - Detector ( $\lambda$ , nm): UV (210) Mobile phase: 50 mM KH <sub>2</sub> PO <sub>4</sub> , H <sub>3</sub> PO <sub>4</sub> (pH 2.8) Type of elution: Isocratic Injection volume (µL): 5 Flow rate (mL/min): 0.6 Column temp: (°C): 20 Run time (min): 25	LOD: - LOQ: -	[58]
Taro corms	Total and soluble oxalates	Column: Rezex ROA ion exclusion organic acid column ( $300 \times 7.8 \text{ mm}$ ; 8 µm particle size) Guard column: Cation H-guard column Detector ( $\lambda$ , nm): DAD (210) Mobile phase: $60\% \text{ H}_2\text{SO}_4$ ( $0.005 \text{ N}$ ) and $40\% \text{ CH}_3\text{CN}$ Type of elution: Isocratic Injection volume ( $\mu$ L): 5 Flow rate (mL/min): 0.5 Column temp: (°C): 50 Run time (min): 20	Linearity range: 6.25 × 10 <sup>-3</sup> –0.2 mg/mL Determination coefficient (r <sup>2</sup> ): 0.999 LOD: - LOQ: -	[68]
Mexican vegetables	Total and soluble oxalates	Column: Sphereclone ODS-column ( $250 \times 4.6 \text{ mm}$ ; 5 µm particle size) Guard column: - Detector ( $\lambda$ , nm): UV- <i>Vis</i> (215) Mobile phase: H <sub>2</sub> SO <sub>4</sub> (1.8 µM) in distilled water (pH 2.6) Type of elution: Isocratic Injection volume (µL): 20 Flow rate (mL/min): 0.4 Column temp: (°C): - Run time (min): -	Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[56]

Matrix	Anal	ytes	Chromatographic Conditions		Validation Data	References
Beans, lentils, sweet potato, and others	Total and solu	Mobile phase: Ad	nPac AS4A anion exchange column (250 size) Guard column: - Detector (λ, nm): Amperometric queous EDTA solution (2.0 g/L) adjuste NaOH Type of elution: Isocratic Injection volume (μL): - Flow rate (mL/min): - Column temp: (°C): - Run time (min): -		Linearity range: - Determination coefficient (r <sup>2</sup> ): - LOD: - LOQ: -	[63]
	ž	r; LOD—limit of detection; LOQ—limit o rophotometric conditions for the ana		"—information not availa	ble.	
Matrix	Analytes	Reaction <sup>(a)</sup>	Spectrophotometric Conditions	Standard Solution	Validation Data	References
Fruits, beverages, and vegetables	Soluble oxalates	Oxidation of oxalate by oxygen in the presence of oxalate oxidase, forming hydrogen peroxide, which will form a monoperoxo complex: $(COOH)_2 + O_2 \rightarrow 2CO_2 + H_2O_2$ $TiO(tpypH_4)^{4+} + H_2O_2 \rightarrow$ $TiO_2(tpypH_4)^{4+} + H_2O.$	Temperature of reaction (°C): 75 pH of reaction: 3 Time of reaction (min): - Oxalate effect: - Wavelength (nm): 450	Sodium oxalate	Linearity range (µg/mL): 0.067–33.5 Determination coefficient (r <sup>2</sup> ): 0.998 LOD (µg/mL): - LOQ (µg/mL): -	[77]
Spinach	Total oxalates	Oxidation of <b>rhodamine B</b> by potassium dichromate in sulfuric acid.	Temperature of reaction (°C): 90 pH of reaction: - Time of reaction (min): 8 Oxalate effect: Catalyst Wavelength (nm): 555	Oxalic acid	Linearity range (µg/mL): 0.06–40 Determination coefficient (r <sup>2</sup> ): - LOD (µg/mL): 0.02 LOQ (µg/mL): -	[38]
Vegetables	Soluble oxalates	Oxidation of <b>bromophenol blue</b> by potassium dichromate in dilute sulfuric acid media.	Temperature of reaction (°C): 60 pH of reaction: - Time of reaction (min): 10 Oxalate effect: Catalyst Wavelength (nm): 600	Oxalic acid	Linearity range (µg/mL): 0.1–8.0 Determination coefficient (r <sup>2</sup> ): 0.998 LOD (µg/mL): 0.04 LOQ (µg/mL): -	[79]

Matrix	Analytes	Reaction <sup>(a)</sup>	Spectrophotometric Conditions	Standard Solution	Validation Data	References
Spinach	Total oxalates	Oxidation of <b>brilliant cresyl blue</b> by potassium dichromate in acidic media.	Temperature of reaction (°C): 80 pH of reaction: - Time of reaction (min): - Oxalate effect: Catalyst Wavelength (nm): 625	Oxalic acid	Linearity range (µg/mL): 0.020–4.70 Determination coefficient (r <sup>2</sup> ): 0.996 LOD (µg/mL): 0.005 LOQ (µg/mL): -	[37]
Spinach	Soluble oxalates	Oxidation of <b>safranine</b> by potassium dichromate in dilute sulfuric acid media.	Temperature of reaction (°C): 60 pH of reaction: - Time of reaction (min): - Oxalate effect: Catalyst Wavelength (nm): 530	Oxalic acid	Linearity range (µg/mL): 0.10–10.0 Determination coefficient (r <sup>2</sup> ): 0.998 LOD (µg/mL): 0.08 LOQ (µg/mL): -	[31]
Spinach	Soluble oxalates	Oxidation of Mn(II) to MnO <sub>4</sub> <sup>-</sup> by potassium periodate: MnSO <sub>4</sub> + KIO <sub>4</sub> $\rightarrow$ MnO <sub>4</sub> <sup>-</sup> + 2IO <sub>3</sub> <sup>-</sup>	Temperature of reaction (°C): 35 pH of reaction: - Time of reaction (min): 18 Oxalate effect: Catalyst Wavelength (nm): 525	Sodium oxalate	Linearity range (µg/mL): 0.05–1.25 and 0.05–1.75 Determination coefficient (r <sup>2</sup> ): 0.998 (for the range of 0.05–1.25 µg/mL) LOD (µg/mL): 0.027 and 0.005 LOQ (µg/mL): -	[80]
Spinach and mushrooms	Total oxalates	Oxidation of <b>pyrocathecol violet</b> with potassium dichromate in acidic media.	Temperature of reaction (°C): 30 pH of reaction: - Time of reaction (min): - Oxalate effect: Catalyst Wavelength (nm): 450	Oxalic acid	Linearity range (μg/mL): 0.08–1.30 Determination coefficient (r <sup>2</sup> ): 0.9993 LOD (μg/mL): 0.07 LOQ (μg/mL): -	[74]
Spinach	Soluble oxalates	Oxidation of <b>Victoria blue B</b> by potassium dichromate in dilute sulfuric acid media.	Temperature of reaction (°C): 60 pH of reaction: - Time of reaction (min): 9 Oxalate effect: Catalyst Wavelength (nm): 610	Oxalic acid	Linearity range ( $\mu$ g/mL): 0.06–9.0 Determination coefficient (r <sup>2</sup> ): $\geq$ 0.996 LOD ( $\mu$ g/mL): 0.12 LOQ ( $\mu$ g/mL): -	[81]
Spinach and mushrooms	Soluble oxalates	Oxidation of iodide by bromate in acidic media catalyzed by iron(II) in the presence of oxalate ion as activator: $BrO_{3-} + 9I^- + 6H^+ \rightarrow 3I_{3-} + 3H_2O$ $+ Br^-$ .	Temperature of reaction (°C): 20 pH of reaction: 5 Time of reaction (min): - Oxalate effect: Activator Wavelength (nm): 352	Sodium oxalate	Linearity range (µg/mL): 0.10–7.0 Determination coefficient (r <sup>2</sup> ): 0.998 LOD (µg/mL): 0.08 LOQ (µg/mL): -	[54]
Spinach, beetroot, and mushrooms	Soluble oxalates	Oxidation of <b>Victoria blue 4R</b> by dichromate in acidic media.	Temperature of reaction (°C): 25 pH of reaction: 4 Time of reaction (min): - Oxalate effect: Catalyst Wavelength (nm): 615	Oxalic acid	Linearity range (µg/mL): 2.0–180 Determination coefficient (r <sup>2</sup> ): 0.995 LOD (µg/mL): 0.7 LOQ (µg/mL): -	[82]

Matrix	Analytes	Reaction <sup>(a)</sup>	Spectrophotometric Conditions	Standard Solution	Validation Data	References
Tap water	Soluble oxalates	Reduction of <b>copper(II) complex</b> to copper (I) by oxalate ion.	Temperature of reaction (°C): RT pH of reaction: - Time of reaction (min): 10 Oxalate effect: - Wavelength (nm): 533	Sodium oxalate	Linearity range (µg/mL): 0.1–2.0 Determination coefficient (r <sup>2</sup> ): - LOD (µg/mL): - LOQ (µg/mL): -	[75]
Spinach and mushrooms	Soluble oxalates	Oxidation of <b>crystal violet</b> by potassium dichromate in sulfuric acid media.	Temperature of reaction (°C): 20 pH of reaction: - Time of reaction (min): - Oxalate effect: Catalyst Wavelength (nm): 630	Sodium oxalate	Linearity range (µg/mL): 0.2–1.8 and 1.8–5.5 Determination coefficient (r <sup>2</sup> ): 0.993 and 0.994 LOD (µg/mL): 0.05 LOQ (µg/mL): -	[78]
Vegetables	Total and soluble oxalates	Ion association of stable anionic complex, which is produced by the reaction of oxalate with Mo(VI) and with Toluidine blue (TBH <sup>2+</sup> ): $(MoO_4)^{2-} + Ox^{2-} \leftrightarrow [MoO_3(Ox)]^{2-}$ $+ H_2O$ $[MoO_3(Ox)]^{2-} + TBH^{2+} \leftrightarrow$ <b>TBH<sup>2+</sup>[MoO_3(Ox)]</b> or $[MoO_3(Ox)]^{2-} + 2TB^+ \leftrightarrow$ <b>(TB)<sub>2</sub>[MoO_3(Ox)]</b> .	Temperature of reaction (°C): - pH of reaction: 6 Time of reaction (min): - Oxalate effect: - Wavelength (nm): 627	Sodium oxalate	Linearity range (µg/mL): 0.0012- 0.012 and 0.0120.240 Determination coefficient (r <sup>2</sup> ): 0.9974 and 0.9915 LOD (µg/mL): 0.00036 LOQ (µg/mL): 0.0012	[30]
Vegetables and aromatic plants	Soluble oxalates	Interaction of oxalates with a sensor material of Si-Ti/ECR, silica-titania xerogel with eriochrome cyanine R.	Temperature of contact (°C): - pH of contact: 3 Time of contact (min): 15 Oxalate effect: - Wavelength (nm): 570	Oxalic acid	Linearity range (µg/mL): 35–900 Determination coefficient (r <sup>2</sup> ): 0.9982 LOD (µg/mL): 10.5 LOQ (µg/mL): 35	[41]
Spinach and mushrooms	Soluble oxalates	Reaction of Reactive blue 4-Cu <sup>2+</sup> with oxalate: <b>RB4-Cu<sup>2+</sup></b> + $(C_2O_4)^{2-}$ $\rightarrow$ CuC <sub>2</sub> O <sub>4</sub> + <b>RB4</b> .	Temperature of reaction (°C): - pH of reaction: 5–7 Time of reaction (min): - Oxalate effect: - Wavelength (nm): 607	Potassium oxalate	Linearity range (µg/mL): 0.29–8.21 Determination coefficient (r <sup>2</sup> ): 0.9983 LOD (µg/mL): 0.10 LOQ (µg/mL): 0.34	[76]
Vegetables and fruits	Total oxalates	Reduction of <b>hexavalent chromium</b> by oxalic acid in presence of Mn(II) as a catalyst.	Temperature of reaction (°C): 25 pH of reaction: 3 Time of reaction (min): 60 Oxalate effect: - Wavelength (nm): 350	Potassium oxalate	Linearity range (μg/mL): 1.660–332.4 Determination coefficient (r <sup>2</sup> ): 0.997 LOD (μg/mL): 0.20 LOQ (μg/mL): 0.66	[58]

<sup>a</sup> Chemical compounds in **bold** are the substances whose absorbance is measured in each reaction. LOD—limit of detection; LOQ—limit of quantification. "-"—information not available.

# 3.2.1. HPLC Conditions

Regarding HPLC conditions for the determination of oxalates (Table 3), ion exchange, ion exclusion, and reversed-phase columns are used. Considering ion exchange columns, one of the most referenced columns is IonPac AS4A [1,59,60,62–64,72]. Other authors have used different ones: propylamine anion exchange column (25 cm  $\times$  4.6 mm; 6  $\mu$ m particle size), diethylaminoethyl (DEAE) anion exchange column ( $250 \times 4.6$  mm; 5 µm particle size), and Hamilton HA-X8.00 column with strong anion exchange resin in the sulfate form  $(255 \times 5 \text{ mm}; 7-10 \text{ }\mu\text{m} \text{ particle size})$  [84]; Alltech All-Sep anion exchange column (100  $\times$  4.6 mm; 7  $\mu$ m particle size) [17] and Shodex IC SI-90 anion exchange column ( $250 \times 4$  mm; 9 µm particle size) filled with KanK-ASt ( $120 \times 5$  mm, 14 µm particle size) [71]. For ion exclusion, it seems that an ion exclusion column ( $300 \times 7.8 \text{ mm}$ ) is frequently the selected one, either from Bio-Rad [32,34,40,61,70] or Rezex [65–68] brands. Regarding reversed-phase columns, the frequent use of  $250 \times 4.6$  mm columns with different particle sizes has been reported [35,36,39,55,56,58,73,84], with 5 µm being the most common. Additionally, a radial compression column with C-18 or C-8 functionality (100 mm; 10 µm particle size) [84] has also been used. Some consumables listed in Table 3 may no longer be available in the market or have been replaced by columns with enhanced characteristics that contribute to the improvement of the analysis.

These columns are chosen considering different separation modes of HPLC. Ion chromatography is an effective method for the determination of oxalate ions because oxalic acid is a strong acid, giving away its protons and becoming negatively charged [71].

When using ion exclusion chromatography, the dissociated functional groups of the ion exchange resin present in the stationary phase have the same charge signal as the oxalate ion and, thus, it is repulsed and eluted [85]. In contrast, ion exchange chromatography is based on the exchange of oxalate ions with the counter-ions, which are anions of the ionic groups attached to the solid support being strongly retained [29]. Reversed-phase chromatography uses a non-polar stationary phase (e.g., C-18) with a polar mobile phase, being the most popular HPLC mode. The separation mechanism is based on polarity and hydrophobic/hydrophilic interactions between oxalates and these two phases. Furthermore, the use of a guard column, a smaller column applied before the analytical column to protect it from impurities present in samples and enhancing the lifetime of the main column, is commonly observed [29]. For the determination of oxalates, C-18 [55], amino [35], IonPac AG4A [59,60], and cation-H [32,34,40,65–68] guard columns are mentioned.

The main used type of elution for oxalate determination is undoubtedly isocratic. In a study using different elution programs for the separation of various organic acids, including oxalic acid, a gradient of  $(NH_4)_2SO_4$  and  $MgSO_4$  resulted in a rise in the baseline of the chromatograms [84]. Thus, using an isocratic elution promotes better results.

For the measurement of oxalates, there is a wide variety of mobile phases. For ion exchange columns, a carbonate and sodium bicarbonate solution for a conductivity detector [17,59,60,71] and an aqueous EDTA solution for amperometric detection are frequently utilized [1,63,64,72]. Some other mobile phases have been reported like 0.15 M  $NaH_2PO_4$  (pH = 4.2); 0.30 M  $NH_4H_2PO_4$  adjusted to pH 6.5 with concentrated ammonium hydroxide; 0.5 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> adjusted to pH 7.25 with ammonium hydroxide; 0.3 M  $(NH_4)_2SO_4$  containing 10% methanol and 0.1–1.5 M MgSO<sub>4</sub> solutions for the gradient analysis [84]. When using ion exclusion columns, sulfuric acid is always the chosen mobile phase. However, Zulkhairi et al. [68] used a mixture of this acid with acetonitrile. For reversed-phase columns, dihydrogen phosphate ( $H_2PO_4^-$ ) is frequently used. This ion can be used alone [35,84] or combined in a mixture with tertiary butyl alcohol (TBA) [55] or tetrabutylammonium hydrogen sulfate [36,39] and it is frequently buffered with phosphoric acid (*circa* pH 2) [36,55,58,84]. Other authors used different solutions as mobile phase, for example, tetrabutylammonium chloride [73] and sulfuric acid [56]. In addition, buffer solutions have a big impact on the retention of analytes [29], so there is frequent use of these solutions to maintain a stable pH. For reversed-phase and ion exclusion columns, acidic buffers from approximately pH 2–3 are commonly used. In contrast, for ion exchange columns, higher pH values are allowed (*circa* pH 5).

The flow rate varies from 0.4–2 mL/min, with 1 mL/min or less being more common. When it is mentioned, column temperature ranges from room temperature to 80 °C. However, usually lower temperatures than 80 °C, such as 20 °C [58], 25 °C [65,66], 30 °C [73], 33 °C [71], or 50 °C, are used [58,68,70]. The most-reported injection volumes are 5 or 20  $\mu$ L.

The most usual detectors are ultraviolet (UV) or ultraviolet–visible (UV–Vis), with 210 nm as the most common wavelength. In contrast, some authors report measurements at 206 [35,84], 215 [40,56], and 220 nm [39,55]. The choice of wavelength is usually made considering the maximum absorbance of analytes and, thus, maximum sensitivity. However, there are other parameters which are also taken into account for this choice, such as analysis time. Other used detectors are diode array detectors (DADs) or photodiode array (PDA) [68], which also detect absorption in the UV to Vis region but can scan the entire range [86]; refractive index (RI) detectors [84]; conductivity detectors [17,59,60,71]; and amperometric detectors [1,63,64,72]. These instruments are chosen depending on the type of HPLC analysis, considering its analytes. For example, authors who employed ion exchange chromatography for oxalate measurement have utilized conductivity detectors, whereas researchers who applied the HPLC–enzyme reactor method (HPLC-ER) used an amperometric detector. This last method is based on the chromatographic separation of oxalate combined with enzymatic conversion to hydrogen peroxide by oxalate oxidase and its amperometric detection [72,87].

#### 3.2.2. Spectrophotometric Conditions

Contrasting with HPLC, studies that describe spectrophotometric methods analyze much fewer samples in quantity and variety, whereas authors who applied HPLC present a large number of results of various foods. In the measurement of oxalates by spectrophotometry (Table 4), spinach and mushrooms are frequently studied. These observations might be due to HPLC's ability for automatization and spectrophotometry's laborious procedures when applied to a large number of samples.

Spectrophotometric conditions are quite peculiar. The studies analyzed for this matter include catalytic/kinetic methods. In other words, these procedures follow a spectrophotometric reaction which has some kind of oxalate intervention, either as a reagent [30,41,58,75–77], catalyst [31,37,38,74,78–82], or activator [54]. Some manuscripts mention other studies in which oxalate acts as an inhibitor [54,78]. These are indirect methods since they do not measure the absorbance of oxalate, but by measuring the absorbance of a substance within a system where oxalate has influence, it is possible to extrapolate oxalate content.

Most of these systems are redox reactions, using a wide range of different reagents (Table 4). For oxalate measurement, the most common type of redox reaction is oxidation in acidic media. However, there are a few reductions reported [58,75]. The use of potassium dichromate as an oxidant agent is frequently observed, but oxygen [77] and bromate [54] have also been reported. In contrast, oxalate has been used as a reducing agent [58,75].

In addition, other more complex procedures exist. For example, Matsubara et al. [77] studied a spectrophotometric method based on the oxidation of oxalate by oxygen in the presence of oxalate oxidase, forming hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) which forms a complex with TiO(tpypH<sub>4</sub>)<sup>4+</sup> and TiO<sub>2</sub>(tpypH<sub>4</sub>)<sup>4+</sup> that absorbs at 450 nm. Furthermore, Mo(VI) can form a stable complex with oxalate,  $[MoO_3(Ox)]^{2-}$ , which can have different forms in a solution in the pH range of 2–7,  $[Mo_2O_5(Ox)_2]^{2-}$ , or  $[Mo_2O_5(OH)(Ox)_2]^{3-}$ . Subsequently, this anionic complex associates with Toluidine blue (TBH<sup>2+</sup>) which has a maximum peak absorbance at 627 nm [30]. A new sensor material for solid-phase spectrophotometric determination of food oxalates was also developed. This method is based on the interaction of oxalates with a material of silica–titania xerogel with eriochrome cyanine R (Si-Ti/ECR) which causes sensor material discoloration, with absorbance being used as an analytical signal [41]. According to Tavallali et al. [76], the reaction of Reactive blue 4 (RB4)-Cu<sup>2+</sup> with

oxalate can be used for oxalate determination in food samples. The addition of oxalate to the RB4-Cu<sup>2+</sup> complex increased the absorption band intensity at 607 nm and changed the color from sky blue to dark blue due to the regeneration of RB4 by the chelation of oxalate with Cu<sup>2+</sup>, since the binding constant of Cu<sup>2+</sup> with oxalate is larger than that of Cu<sup>2+</sup> with RB4.

Temperature, pH, and time of reactions are highly specific to each reaction, being parameters optimized before selecting the final procedure. The values of these parameters are chosen considering sensitivity and reproducibility [31,78]. Temperature, pH, and time of these reactions can range from approximately 20–90 °C, 3–7, and 8–60 min, respectively.

Reactions are monitored by measuring the absorbance of the reagents or products which are chromophore substances, such as crystal violet, Victoria blue, and brilliant cresyl blue, at maximum wavelength,  $\lambda_{max}$ , the wavelength whose absorbance is maximum and producing maximum sensitivity. This measurement is always in the UV–*Vis* region (200–800 nm), mainly in the visible region (400–800 nm), since analyzed compounds are frequently colored and absorb this kind of light. For the construction of calibration graphs, oxalic acid as the standard solution was mainly used, but sodium oxalate [30,54,75,77,78,80] and potassium oxalate [58,76] were also reported.

By reading spectrophotometric methods for the determination of oxalate content, it can be concluded that validation is a frequent concern. Thus, linearity range, determination coefficient ( $r^2$ ), and limit of detection (LOD) are often presented parameters (Table 4).

In Figure 2, a summary of the most commonly used conditions for extraction for HPLC and spectrophotometry analysis is provided. However, according to the literature, a wide variety of techniques and conditions are applied depending on the type of food, which makes it difficult to define a single method to measure oxalate in foods.



**Figure 2.** Summary of the most commonly used conditions for extraction for HPLC and spectrophotometry measurement of oxalate in foods. <sup>(a)</sup>—These parameters are highly specific to each reaction (Table 4).

# 4. Oxalate Occurrence in Foods

Oxalate contents (mg/100 g FW) of various foods measured by chromatographic and spectrophotometric methods, such as in fruits, vegetables, mushrooms, legumes, pseudoce-reals, and aromatic plants, are collected in Table 5 in alphabetical order. Other studies were

also taken into consideration; however, their oxalate contents were not included in this table because they were presented in dry weight and it was not possible to convert them into fresh weight owing to the lack of moisture or dry matter values [57,58,60,62,66,68].

Table 5. Occurrence of oxalate (mg/100 g of fresh weight) in foods.

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References
Amaranth Amaranthus mangostamus	-	Raw	1510.8	835.1	675.7	[40]
Apple Malus domestica	1	Granny Smith, raw	3.5	1.8	1.7	
Apple Malue culturetrie	26	- Cox Kent, raw	_	0.4	0.6	-
Apple Malus sylvestris	28	- Cox Kent, raw	1	-	0.6	[1]
Apricot Prunus armeniaca	2	Raw	6.8	1.9	4.9	
Artichoke <i>Cynara</i> scolymus	1	Boiled	6.8	6.8	0	
Arugula Eruca vesicaria	10	Raw	_	ND	-	[59]
Asparagus Asparagus officinalis	3	Boiled	2.6	0.9	1.7	[1]
Asparagus chicory Cichorium intybus	10	Raw	-	1	-	[59]
Aubergine <i>Solanum</i>	3	Green, long, and raw; boiled	55; 38	45; 19	10; 19	[61]
	-	Raw	54.4	53.7	0.7	[40]
melogena (or eggplant)	1	Boiled	12.8	4.8	8	
-	2	Raw	16.2	15.7	0.5	[1]
Avocado Persea gratissima	2	Raw	1.3	1.3	0	
Bamboo shoot Bambusa	3	Cultivated and raw; boiled	222; 93	163; 51	60; 42	[(1]
spp.	3	Pickled and raw; boiled	71; 51	23; 10	47; 41	[61]
Banana Musa paradisiaca	7	Raw	6.8	0.7	6.1	[1]
	2	Preserved white	54.2	1.9	52.3	- [1]
-	2	White, seeds, dry	547.9	38.8	509.1	[63]
-	3	White, raw; boiled	158; 47	52; 12	106; 34	[26]
	3	Red, raw; boiled	113; 72	37; 22	76; 50	[36]
	1	Quail, seeds, dry	176.7	16.9	159.8	[63]
Bean Phaseolus vulgaris	3	Red kidney, seeds, and raw; boiled	91; 32	26; 10	65; 27	[61]
	1	Kidney	-	32	-	[79]
-	1	Red kidney	13.9	1.5	12.4	[1]
	1	Red kidney, seeds, dry	74.6	4.8	69.8	[63]
-	1	Green, raw	65.2	8.4	56.8	

	Tał	<b>ble 5.</b> Cont.					
Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References	
	3	Raw; boiled	67; 52	45; 38	22; 15	[36]	
-	5	Raw	-	1431	-	[82]	
Beetroot Beta vulgaris	2	Boiled	36.9	16.3	20.6	[1]	
-	6	Raw	-	74,9	-	[59]	
	-	Raw; boiled	45.6; 76.0	38.6; 72.3	7; 3.7	[32]	
Bellflower root Platycodon grandiflorum	-	Raw	ND	ND	-	[40]	
Bilberry Vaccinium myrtillus	3	Raw	1.5	0.1	1.4		
Blackcurrent Ribes nigrum	1	Raw	19	3	16	[1]	
Bramble (blackberry) Rubus fruticosus	4	Raw	29.2	0.9	28.3	-	
	-	Raw	ND	ND	-	[40]	
Broccoli Brassica oleracea	2	Boiled	1.4	1.1	0.3	[1]	
Broccoll Brassica oleracea –	5	Raw	-	0,5	-	[59]	
	-	Raw; boiled	16.1; 10.1	11.6; 6.6	4.5; 3.5	[32]	
Broccoli raab Brassica rapa	10	Raw	-	0,1	-	[59]	

Broccoli Brassica oleracea	2	Boiled	1.4	1.1	0.3	[1]
Broccon Brussicu oleruceu	5	Raw	-	0,5	-	[59]
	-	Raw; boiled	16.1; 10.1	11.6; 6.6	4.5; 3.5	[32]
Broccoli raab Brassica rapa	10	Raw	-	0,1	-	[59]
Brussel sprout Brassica oleracea	3	Boiled	1.2	0.8	0.4	[1]
Burdock Arctium lappa	1	Raw	64.8	62.7	2.1	[40]
	-	Raw	ND	ND	-	[40]
Cabbage Brassica oleracea	3	Raw; boiled	7;5	ND; 4	7; <dl< td=""><td>[61]</td></dl<>	[61]
	3	Raw	-	ND	-	[59]
	3	Raw; boiled	29; 12	24; 7	5; 5	[61]
	3	Raw; boiled	49; 26	28; 14	21; 12	[36]
	-	Raw	16.4	16.2	0.2	[40]
Carrot Daucus carota	24	Raw	17.8	9.0	8.8	[1]
	1	Boiled	4.9	2.3	2.6	[1]
	7	Raw	-	12	-	[59]
	-	Raw; boiled	35.6; 32.3	22.6; 19.3	13.0; 13.0	[32]
	1	D	-	0.3	0.1	[1]
Cauliflower Brassica	2	Raw	0.4	-	0.1	[1]
oleracea	4	Raw	-	ND	-	[59]
	3	Raw; boiled	27; 8	ND; ND	27; 6	[61]
	-	Raw	23.2	ND	23.2	[40]
Celery Apium graveolens	1	Canned	6.7	3.5	3.2	[1]
	6	Raw	-	0,5	-	[59]

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	Reference
Chanterelles Cantharellus cibarius	2	Canned	0.5	0.5	0	[1]
Cherry Prunus avium	5	Sweet, raw	2.4	1.3	1.1	
Chickpeas Cicer arietinum	2	Seeds, dry	14.3	13.7	0.6	[63]
Chinese convolvulus Lpomoea reptans	3	Raw; boiled	156; 135	21; 7	135; 128	
Chinese kale Brassica oleracea	3	Raw; boiled	23; 7	ND; ND	22; 7	[61]
Chinese cucumber Momordica charantia	3	Raw; boiled	71; 56	57; 22	14; 34	
Coriander Coriandrum	-	Dried	40.6	ND	40.6	[33]
sativum	4	Raw, air dried	40.5	ND	40.5	[34]
Cowpea Vigna unguiculata	3	Seeds and raw; boiled	ND; 5	ND; ND	ND; 5	[61]
Cress Nasturtium officinale	1	-	ND	ND	-	[1]
Crown daisy Chrysanthemum coronarium	-	Raw	96.0	58.8	37.2	[40]
Cucumber Cucumis	-	Raw	ND	ND	-	
sativus	1	Raw	0.4	0.3	0.1	
Elderberry Sambucus nigra	4	Black, raw	72.1	7.1	65	[1]
Endive Cichorium endivia	5	Raw	-	0.2	-	[59]
Fava beans Vicia faba	1	Seeds, raw	1.3	0.9	0.4	[63]
	2	Boiled	5.3	3.3	2	[1]
Fennel Foeniculum vulgare	1	Raw	19.7	17.2	2.5	[1]
0	11	Raw	-	12.4	-	[59]
_	3	Raw	20.5	3.3	17.2	
Fig Ficus carica	2	Dried	-	5.5	89.6	[1]
	1	Diled	95.1	-	09.0	
Garlic Allium sativum	-	Raw	ND	ND	-	[40]
	4	Raw	_	ND	_	[59]
Green amaranth	-	Leaves, dried	1939.9	901.7	1038.1	[33]
Amaranthus viridis	4	Raw, air dried	1940.8	902.3	1038.7	[34]
Green peas Pisum sativum	2	Seeds, dry	3.3	2.7	0.6	[63]
Green onion Allium cepa	3	Raw	-	ND	-	[59]
Green pepper Capsicum annuum	-	Raw	31.0	27.5	3.5	[40]

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References	
Gooseberry Grossularia	1	Red	21.6	3.2	18.4		
uva crispa	1	Green, raw	27	3.1	23.9	-	
Granadilla Passiflora edulis	1	Raw	1	0.6	0.4	[1]	
Grape Vitis vinifera	9	Green, raw	1.7	0.6	1.1	-	
Huauzontle Chenopodium Nuttalliae	3	Raw; boiled	162.0; 97.7	135.0; 97.7	27.1;0	[56]	
Ivy gourd Coccinia grandis	3	Raw; boiled	36; 24	10; 5	29; 19	[61]	
Kale Brassica oleracea	-	Raw	ND	ND	-	[40]	
Kiwi fruit Actinida	-	-	-	4.50	-	[77]	
chinesis	6	Raw	23	2.4	20.6	[1]	
Killin in the start	1	Boiled	0.7	0.7	0	- [1]	
Kohlrabi Brassica oleracea	3	Raw	-	ND	-	[59]	
Leaf chicory Cichorium	-	Raw	47.6	42.3	5.3	[40]	
intybus	8	Raw	-	ND	-	[59]	
Leek Allium porrum	2	Raw	17.0	9.4	7.6	[1]	
Leek Allium tuberosum	-	Raw	48.6	45.1	3.5	[40]	
Lemon Citrus medica	2	Raw	3.1	0.5	2.6	[1]	
Lentil Lens culinaris	2	Brown, seeds, dry	24.0	14.3	9.7	[(0]	
Lentil Lens culturis	2	Red, seeds, dry	13.8	7.8	6	- [63]	
Lentil Lens esculenta	1	Dried	13.3	1.9	11.4	[1]	
	-	Raw	40	ND	40	[40]	
Lettuce <i>Lactuca sativa</i>	2	Raw	0.3	0.3	0	[1]	
-	16	Raw	-	ND	-	[59]	
Lettuce Valeriana locusta	1	Raw	1.3	0.9	0.4	[1]	
Licorice Glycyrrhiza glabra	2	Root	3569.3	165	3404.3	[63]	
Lime Citrus auranttifolia	1	Raw	7.5	0.5	7	[1]	
Lotus root Nelumbo nucifera	-	Raw	ND	ND	-	[40]	
Mandarin Citrus nobilis	2	Raw	8.5	0.3	8.2		
	4	P	-	0.5	1.1	-	
Mango Magnifera indica	5	Raw	1.6	-		[1]	
Mirabelle Prunus domestica syriaca	4	Raw	8.1	0.2	7.9		
Mung bean Vigna radiata	3	Seeds and raw; boiled	24; 5	12; ND	12; 4	[61]	

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References
	5	Raw	-	320	-	[78]
-	1	Raw	-	41	-	[79]
-	5	Raw	-	36	-	[82]
Mushroom Agaricus	1	Raw		-	-	[74]
bisporus	5	Raw	-	326	-	[54]
-	3	Canned	0.7	0.4	0.3	
-	1	Boiled	0.5	0.1	0.4	- [1]
	-	Raw	-	482	-	[76]
Muskmelon Cucumis	6		-	0.9		
melo	7	Raw	1	-	- 0.1	[1]
New Zealand Spinach Tetragonia tetragonioides	-	Raw; boiled	1764.7; 1322.6	364.6; 129.3	1400.1; 1193.3	[32]
Okra Abdelmoschus esculentus	1	Raw	317.2	56.3	260.9	[63]
01: 01	2	Green, canned	45.7	1.2	44.5	[4]
Olive Olea europaea	3	Black, canned	13.9	1.6	12.3	- [1]
	-	Raw	ND	ND	-	[40]
Onion Allium cepa	2	Raw	1.7	1.6	0.1	[1]
	6	Raw	-	ND	-	[59]
	-	Dried	29.1	ND	29.1	[33]
Onion stalks <i>Allium cepa</i>	4	Raw, air dried	29.3	ND	29.3	[34]
Orange Citrus sinsensis	1	Raw	1.8	0.2	1.6	
Oval kumquat Fortunella margarita	1	Raw	3.5	0.8	2.7	[1]
Papaloquelite Porophyllum ruderale	3	Raw	15.8	15.8	ND	[56]
Demons Cavias namena	3	Raw; boiled	5; 11	ND; ND	ND; 8	[61]
Papaya Carica papaya	3	Raw	1.3	0.5	0.8	[1]
Parsley	3	Raw	-	782	-	[41]
Parsley Petroselinum	1		-	76	<i>(</i> 2	[1]
sativum	2	-	136	-	- 60	[1]
Parsley Petroselinum	-	Raw	270.7	72.0	198.8	[40]
crispum	8	Raw	-	0.5	-	[59]
	1	Canned	6.2	6.2	0	
Pea Pisum sativum	1	Boiled	0.2	0.2	0	-
-	2	Green, dried	ND	ND	-	-
Peach Prunus persica	3	Raw	2.5	0.2	2.3	-
	10	_	-	0.9		- [1]
	11	Raw	2.7	-	- 1.8	
Pear Pyrus communis	2		-	3.8		-
	3	Peeled	3.7			

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References
Persimmon Diospyros kaki	-	-	-	2.61	-	[77]
Pineapple Ananas	2	Preserved without	-	0.9	4	
comosus	3	sugar	4.9	-	- 4	
Plum Prunus domestica	8	Raw	1.7	0.5	1.2	-
	3	Boiled	24.3	12.8	11.5	[1]
-	1	Baked	13.0	11.7	1.3	- [1]
- Potato <i>Solanum tuberosum -</i> - -	1	Deep fried	26.9	17.0	9.9	-
	1	Chips	47.0	45.8	1.2	-
	2	Raw	17.1	13.0	4.1	-
	3	Raw	-	0,4	-	[59]
	-	Raw	ND	ND	-	[40]
Pumpkin Cucurbita pepo	1	Raw	ND	ND	-	[1]
Purple amaranth	-	Leaves, dried	1355.3	594.9	760.4	[33]
Amaranthus cruentus	4	Raw, air dried	1353.7	594.2	759.5	[34]
Quelite Chenopodium album	3	Raw; boiled	110.0; 59.8	72.4; 31.0	37.6; 28.9	[56]
Radicchio Cichorium intybus	3	Raw	-	ND	-	[59]
	-	White	-	1.46	-	[77]
Radish Raphanus sativus	1	Raw, red	1.7	1.4	0.3	[4]
-	1	Raw, white	ND	ND	-	- [1]
	-	Dried	12.3	ND	12.3	[33]
- Radish leaves	4	Raw, air dried	12.3	ND	12.3	[34]
-	6	Raw	-	ND	-	[59]
~	-	Raw	ND	ND	-	[40]
Radish roots -	6	Raw	-	ND	-	[59]
Raspberry Rubus idaeus	4	Raw	18.9	3.4	-	
Redcurrant Ribes rubrum	4	Raw	19.8	4.9	-	
Rhubarb Rheum	1	_	-	380		- [1]
rhabarbarum	2	Raw	1235	-	855	
Rhubarb petioles	4	Raw	1080	-	-	[55]
Rhubarb stalks	-	Raw; boiled	986.7; 756.3	287.3; 80.7	699.4; 675.6	[32]
Romerito Suaeda torreyana	3	Raw; boiled	94.2; 24.0	67.0; 11.2	27.2; 12.8	[56]
Sauerkraut Brassica oleracea	1	Raw	7.1	7.1	0	[1]
Savoy cabbage Brassica	1	Boiled	3.5	1.3	2.2	
oleracea	2	Raw	_	0.1	-	[59]

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	Reference
Sorrel Rumex acetosa	3	Raw	-	1079	-	[41]
Soybean <i>Glycine max</i>	3	Seeds, raw	204	58	145	[61]
	3	Raw; boiled	497; 224	155; 64	343; 162	[36]
	1	Flakes	218.4	32.6	185.8	[63]
	1	Flour	124	29.9	94,1	
	2	Seeds, dry	276.8	37.9	238.9	-
Soybean sprout	1	Raw	26.5	12.7	13.8	[40]
· •	-	Raw	2350	2000	350	[39]
	-	Raw	460	-	-	[38]
	10	Raw	390	-	-	[37]
Spinach Spinacia oleracea	2	Raw	-	731	-	[80]
	5	Raw	-	762	-	[54]
	3	Raw	-	2.61	-	[81]
	5	Raw	-	415	-	[82]
	2	Raw		-	-	[74]
	3	Raw	-		-	[79]
	-	Raw	-	421	-	[31]
	5	Raw	-	1330	-	[78]
	-	-	-	515	-	[77]
	182	Frozen, cultivated in winter	-	1092.9	-	[73]
		Frozen, cultivated in spring	-	890.3	-	
		Frozen, cultivated in summer	-	752.5	-	
		Frozen, cultivated in fall	-	614.9	-	
	3	Raw; boiled	978; 477	543; 184	435; 293	[36]
	2	Boiled	364	101	263	[1]
	1	Boiled with cream	412	123	289	
	-	Raw	1272.2	1176.1	96.0	[40]
	12	Raw	-	542.6	-	[59]
	-	Raw; boiled	329.6; 154.8	266.2; 90.9	63.4; 63.9	[32]
	1	Raw	1370	90	1280	[35]
	-	Dried	768.4	727.1	41.3	[33]
	4	Raw, air dried	691.7	643.5	48.2	[34]
	-	Raw	-	631	-	[76]
Spring onion Allium fistulosum	-	Raw	33.3	28.2	5.1	[40]

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References
Star fruit Averrhoa carambola	4	Raw	295.4	138.9	156.5	[1]
	-	Raw	160	90	70	[35]
Strawberry Fragaria	8	Raw	2.9	0.9	2	[1]
Sultana Vitis vinifera	1	Dried	8.5	3.2	5.3	[1]
Swamp morning glory Ipomoea aquatica	3	White stems and raw; boiled	79; 56	58; 40	21; 16	[61]
	3	Red stems and raw; boiled	94; 59	61; 23	33; 36	
Sweet potato Ipomea batatas	2	Raw	495.6	76.7	418.9	[63]
	3	Raw	-	ND	-	[59]
Swiss chard Beta vulgaris	6	Raw	874	327	547	[1]
	-	Raw	1458.1	1082.7	375.4	[40]
-	12	Raw	-	207,7	-	[59]
Swiss chard stems	-	Raw; boiled	127.4; 148.4	24.0; 19.4	103.4; 129.0	[32]
Swiss chard leaves	-	Raw; boiled	525.5; 291.1	252.3; 117.7	273.2; 173.4	[32]
	4	Maori cultivar, raw	443.7	204.1	239.5	[34]
Taro leaves Colocasia esculenta	4	Maori cultivar, baked	721.9	367.2	354.7	
	4	Maori cultivar, baked with milk	397.5	173.3	224.2	
	4	Japanese cultivar, raw	424.7	267.0	157.7	
	4	Japanese cultivar, baked	533.9	352.6	181.3	
	4	Japanese cultivar, baked with milk	300.2	144.6	155.6	
Tassel hyacinth Leopoldia comosa	3	Raw	-	ND	-	[59]
Tamarillo <i>Tamarindus</i> indica	1	Raw	19.9	3.7	16.2	[1]
Tomato Lycopersicum esculentum	3	Raw	11	7	4	[61]
	3	Raw	8.5	3.6	4.9	- [1]
	1	Canned, peeled	12.7	3.1	9.6	
Verdolaga Portulaca oleracea	3	Raw; boiled	91.5; 119.8	62.7; 50.1	28.7; 69.7	[56]
Viper's grass Scorzonera hispanica	1	Canned	9.1	6.5	2.6	[1]
Water dropwort Oenanthe javanica	-	Raw	93.0	64.7	28.3	[40]
Watermelon Citrullus lanatus	1	Raw	0.3	0.3	0	[1]

Food	n	Type of Sample	Total Oxalates (mg/100 g FW)	Soluble Oxalates (mg/100 g FW)	Insoluble Oxalates (mg/100 g FW)	References
Winged bean Psophocarpus tetragonolobus	3	Pods and raw; boiled	7;5	ND; ND	4; 5	[61]
Yard long bean Vigna sesquipedalis	3	Green and raw; boiled	38; 29	9; 3	29; 23	
Yellow plum <i>Prunus</i> domestica ssp syriaca	2	Raw	1.4	0.4	1	[1]

ND—not detected; <DL—below limit of detection; n—number of samples; FW—fresh weight.

In general, fruits are considered low-oxalate foods (<30 mg total oxalate/100 g FW), except for star fruit (160 [35] and 295.4 [1] mg total oxalate/100 g FW), elderberry (72.1 mg total oxalate/100 g FW), and dried fig (95.1 mg total oxalate/100 g FW) [1]. However, it is to be noted that for elderberry and dried fig, the majority of oxalate content is insoluble (72.1 and 89.6 mg/100 g FW, respectively) and, thus, less harmful [1].

Regarding vegetables, there is much more variety in the values of oxalate content, ranging from not detected to high amounts of oxalate. For example, raw New Zealand spinach and typical spinach have been reported as extremely high in oxalates with values of 1764.7 [32] and 329.6–2350 mg total oxalate/100 g FW [32–39], respectively. Soluble oxalate concentration in *Spinacia oleracea* has been studied during the cultivation season. Oxalate content was higher in winter (1092.9 mg soluble oxalate/100 g FW) and lowest in fall (614.9 mg soluble oxalate/100 g FW), indicating that higher oxalate content is related to a longer growing period since this compound is an end product of some metabolisms and it increases in the vacuole with plant aging [73]. Also, rhubarb contained high oxalate content (1235 mg total oxalate/100 g FW) [1] as well as Swiss chard with 874 [1] and 1458.1 [40] mg total oxalate/100 g FW, having more in leaves than in stems [32], and sorrel with 1079 mg soluble oxalate/100 g FW [41]. Taro leaves yielded 300.2–721.9 mg total oxalate/100 g FW, depending on the type of cultivar and processing techniques: baking increased oxalate content compared to raw samples, due to concentrating effects, whereas baking with milk (a source of calcium) decreased oxalate content, especially soluble [34]. This decrease happens because oxalate ions can bind with calcium, precipitating and reducing soluble oxalate [45]. Conversely, there were some vegetables with small amounts of oxalate or it was even undetected, e.g., cabbage, broccoli, cauliflower, cucumber, kale, and pumpkin (Table 5).

Regarding legumes, soybean is considered high in oxalates, ranging between 124 and 497 mg total oxalate/100 g FW, depending on the kind of sample analyzed [36,61,63]. For different types of beans, oxalate values ranged between 13.9 and 547.9 mg total oxalate/100 g FW [1,36,61,63], whereas chickpeas [63] and lentils [1,63] yielded low oxalate content (<24 mg total oxalate/100 g FW) and oxalate in cowpea was not detected [61].

In the pseudocereals category, amaranth is considered a high-oxalate food (1510.8 mg total oxalate/100 g FW) [40] as well as green amaranth (*circa* 1940 mg total oxalate/100 g FW) and purple amaranth (*circa* 1354 mg total oxalate/100 g FW) [33,34].

For aromatic plants, two different species of parsley, water dropwort (also known as water celery), and coriander yielded 136 [1] and 270.7 [40], 93 [40], and *circa* 41 [33,34] mg total oxalate/100 g FW, respectively. Licorice was the highest-oxalate food reported in this review with 3569.3 mg total oxalate/100 g FW [63]. In contrast, arugula, cress, garlic, and green onion were some examples of aromatic plants which did not contain detectable oxalate (Table 5).

Differences between oxalate values for the same food, as observed for beans, lettuce, parsley, mushrooms, or spinach, can vary according to growth, ripeness, climate, region,

soil conditions, and time of harvest. In addition to these conditions, which are harder to control, are sample preparation, which can lead to oxalate generation or its incomplete extraction, and analytical methods with different features that can have an impact on oxalate results [1,6,36,65].

It is important to consider the usual amount of consumption of these foods. For example, some aromatic plants, like parsley, contained high oxalate values but their daily intake is naturally much less than 100 g. Also, the type of consumption has to be taken into account because some foods are not generally eaten in the form in which they were analyzed (e.g., raw beans or mushrooms). It is well known that the same food prepared differently (raw, boiled, baked, fried, etc.) can lead to different oxalate results. The reported values for potato (*Solanum tuberosum*) are a great example to evidence the influence of cooking techniques on oxalate content [1]. In general, boiling has been associated with decreased oxalate content, especially soluble oxalates, due to its leaching and thermal degradation [45,47], as observed in spinach, New Zealand spinach, red and white beans, soybean, and rhubarb (Table 5).

### 5. Health Implications of Oxalates

Oxalate has been a concern for human health due to its antinutritive effects and potential nephrotoxicity for a long time [42,43]. In 1989, a fatality from oxalic acid poisoning was reported. A 53-year-old man, with other conditions, had eaten a sorrel soup with 6–8 g of oxalic acid [88]. A lethal dose of oxalic acid for adults was estimated as 10–15 g, although the ingestion of 4–5 g of oxalate was considered the minimum dose able to cause death [51,88]. As antinutrients, oxalates restrict the bioavailability of some nutrients since they can bind to minerals, like calcium, magnesium, or iron, reducing their absorption and use [3,44].

The sources of oxalates in our body can be exogenous or endogenous (Figure 3). Exogenous oxalate sources are mainly plant foods, like vegetables, grains, legumes, and fruits, among others. When these types of foods are ingested, oxalate is absorbed in multiple parts of our gastrointestinal tract, namely the stomach, small intestine, and large intestine. However, the absorption depends on its availability, among other individual features. Insoluble oxalates are excreted in feces since they are less bioavailable and, therefore, pose a lower health risk. In contrast, soluble oxalates are absorbed through the intestines and colon (5–10% of ingested oxalate, under normal conditions), going into the bloodstream.

Since absorption of oxalates is related to the amount of soluble oxalates, which are more bioavailable, a simultaneous consumption of oxalate with calcium or magnesium can reduce its bioavailability and absorption due to the formation and fecal excretion of insoluble salts, lowering the health risk [1,45,47,52,89]. It has been reported that men with less than 755 mg/day of calcium intake had a higher risk of kidney stone formation, whereas men with a median calcium intake or above had a lower risk [3]. Therefore, dietary calcium intake has been inversely associated with kidney stone formation [3,52,90,91].

Also, it has been observed that intestinal absorption of oxalates in individuals with a history of stone formation was expressly higher than in healthy individuals (9.2% and 6.7%, respectively) [1]. Gastrointestinal health influences oxalate absorption as well, with soluble oxalate being excessively absorbed due to intestinal malfunction.

Despite these facts, oxalate is not typically consumed daily in high concentrations and there are other constituents in foods which have a protective role against kidney stone formation, such as phytate, potassium, calcium, and antioxidant phytochemicals like polyphenols [3]. Also, boiling, steaming, soaking, and processing with calcium sources are some procedures to reduce the content of soluble oxalates, the most harmful oxalates [45,52].

Concerning the endogenous production of oxalates, the liver is the primary source. There are different pathways for oxalate production, including the metabolism of protein (through amino acids, like tyrosine, tryptophan, phenylalanine, and hydroxyproline), ascorbic acid, and precursors of oxalate (such as L-glycerate glycollate and glyoxylate) [92,93].



Glyoxylate is an important intermediary product in several reactions and, for its metabolization into oxalate, enzymes like glycolate oxidase and lactate dehydrogenase are needed.

**Figure 3.** Sources, endogenous pathways, and excretion of oxalates in the human body. HOGA1—4hydroxy-2-oxoglutarate aldolase type 1; GRHPR—glyoxylate reductase/hydroxypyruvate reductase; LDHA—lactate dehydrogenase A; ALDH—aldehyde dehydrogenase; GO—glycolate oxidase; AGT—alanine/glyoxylate aminotransferase; DAO—D-amino acid oxidase.

Free oxalates are delivered to the kidney and can be excreted, increasing urinary oxalates, or can chelate with calcium ions there, resulting in calcium oxalate crystals, which can cause serious health issues like kidney stones, also known as nephrolithiasis or urolithiasis (Figure 3) [3,45–48,91].

This crystallization in the kidney infiltrates vessel walls and can lead to renal tubular obstruction, vascular necrosis, and hemorrhage, which can cause anuria, uremia, electrolyte disturbances, or even rupture and kidney failure [48,51]. Calcium oxalate and its relationship with kidney stone formation have been amply studied, with calcium oxalate being one of the most common types of human kidney stone reported, followed by calcium phosphate [46,51,90,94,95].

Hyperoxaluria is a metabolic disease that leads to excessive urinary oxalate excretion (>40–45 mg/day) [89,96], being an indicator of possible kidney stone formation [91]. The most reliable way to assess daily oxalate intake is through 24 h urine collection; however, there are also food frequency questionnaires whose credibility is debated [89,91].

Hyperoxaluria can be divided into primary hyperoxaluria (PH1) and secondary hyperoxaluria (PH2). PH1 is a group of rare autosomal recessive diseases that negatively affect key enzymes of oxalate metabolism, leading to an overproduction of oxalates in the liver [50]. When renal function declines and excess oxalate exists in the bloodstream, a phenomenon known as systemic oxalosis occurs and calcium oxalate crystals deposit in various organs, tissues, and bones [50,96,97]. Severe damage in the eyes, joints [98], myocardium, skin [99], oral tissues [96], and bone marrow [49] is reported. Oxalate can also be associated with acute kidney injury, a tubular obstruction due to calcium oxalate crystal deposition, and with chronic kidney disease progression, but further studies are necessary [89].

Patients with hyperoxaluria, especially PH1, from a clinical point of view, frequently present severe bone pain, pathological fractures, and bone deformations. This is frequently

associated with the fact that calcium oxalate crystals may deposit within bones, tendons, cartilage, and synovium, causing oxalate arthritis. Then, the calcium oxalate crystals may enter into the synovial fluid, where an inflammatory response will arise, leading to joint effusions and arthralgias [100,101].

PH2 results from increased intestinal absorption of dietary oxalates and can also lead to excessive urinary oxalate [48,51,102]. A high intake of foods rich in oxalate, enteric hyperoxaluria, oxalate-degrading mechanisms, and SLC4 and SLC26 ionic exchangers are linked with PH2. Dietary oxalate plays an important role in PH2, contributing up to 72% of the urinary oxalate excreted [17]. Enteric hyperoxaluria is a form of PH2 that is linked with malabsorption syndromes due to disease or resection of the gastrointestinal tract. In foods, oxalate is usually complexed with calcium, resulting in insoluble oxalate, which is difficult to absorb. Nevertheless, in fat malabsorption conditions, the amount of free oxalates can increase, due to the capacity of free fats to bind calcium. Therefore, PH2 is linked with several conditions that cause fat malabsorption, such as inflammatory bowel disease, celiac disease, short bowel syndrome, and bariatric surgery, among others [100].

The gut microbiome plays an important role since some bacterial species can degrade oxalate to obtain carbon and energy and therefore reduce the concentration of oxalates in blood and urine, minimizing the formation of kidney stones [91–93,103,104].

The gut microbiota is usually similar between individuals; however, it can be affected by the age of individuals, by the diet, and by the use of antibiotics, among other factors. Probiotics are defined as "live microorganisms which, when administered in adequate amounts, confer a health benefit on the host" and are being abundantly used as preventive therapeutic agents for several diseases, since they have been implicated in the stabilization of gut microbiota and enhancement of immune responses [105].

The best-known oxalate-degrading microorganism is *Oxalobacter formigenes*, but others are also able to degrade oxalate into carbon dioxide and formate, namely *Escherichia coli, Bifidobacterium* spp., and *Lactobacillus* spp. [93,104]. *O. formigenes* is a Gram-negative anaerobic bacterium isolated from human feces and other animals that utilizes intestinal oxalate as a carbon source, through formyl-CoA transferase and oxalyl-CoA decarboxylase enzymes, and metabolizes the oxalate into carbon dioxide and formate, contributing to regulating oxalate homeostasis [103]. However, its application as a probiotic is limited due to its nutritious requirements, but also because it has less colonization ability and is sensitive to the use of certain antibiotics and drugs [106]. Moreover, the therapeutic use of *O. formigenes* can be compromised, for example, in patients with PH1 and patients with cystic fibrosis. To date, the best conditions (pH, sugar concentration), as well as the adequate amount of these supplements, are not clear and more research is still needed [93].

### 6. Conclusions

Various methods have been employed for the determination of oxalate in foods. Particularly, extraction and analytical conditions of HPLC and UV–*Vis* spectrophotometry were reviewed. Despite having different features, both methods have similar extraction procedures. Among other extraction parameters, temperature remains a controversial question because it can lead to oxalate formation from precursors or failure to dissolve all pre-existing insoluble oxalates. Furthermore, a considerable quantity of different HPLC and spectrophotometry methods were gathered and analyzed, concluding that there is a huge variety of procedures.

This review also compared the oxalate content (mg/100 g FW) of a wide range of foods, measured by HPLC and spectrophotometry. The results showed that spinach, New Zealand spinach, rhubarb, Swiss chard, taro leaves, sorrel, soybean, amaranth, parsley, and licorice contained high oxalate levels and can be considered high-oxalate foods, especially some green leafy vegetables. In contrast, others can be referred to as low-oxalate foods: fruits (except for star fruit, elderberry, and dried fig), cabbage, broccoli, cauliflower, cucumber, kale, pumpkin, chickpeas, lentils, cowpea, arugula, cress, garlic, and green onion.

Nevertheless, there are some procedures to reduce oxalate, in particular soluble oxalate, such as boiling, steaming, soaking, and processing with calcium sources.

Despite a clear relationship between dietary oxalate, calcium oxalate, and kidney stone risk, the connection might be more complex than previously thought due to the impact of cooking techniques, calcium intake, endogenously produced oxalate, and gastrointestinal health. Foods which contain oxalates, such as fruits and vegetables, have a wide range of beneficial compounds that might outweigh possible negative implications on human health. Additionally, systemic oxalosis does not seem to be related to dietary oxalate, but to previous pathologic conditions of individuals such as primary hyperoxaluria. Hence, regular consumption of high-oxalate foods by healthy individuals as a part of a balanced and diversified diet does not appear to cause health issues if daily consumption is from 50–200 mg/day, whereas for individuals susceptible to kidney stone formation dietary modification is crucial for its prevention. For these individuals, it is recommended to limit the consumption of high-oxalate foods to less than 40–50 mg oxalate/day since they can present a health threat in these cases.

**Author Contributions:** Conceptualization, N.S. and T.G.A.; methodology, N.S. and M.A.S.; data collection, N.S. and T.G.A.; writing—original draft preparation, N.S., T.G.A. and M.A.S.; writing—review and editing, M.E.F. and H.S.C.; supervision, M.E.F. and H.S.C. All authors have read and agreed to the published version of the manuscript.

**Funding:** This research received financial support from Foundation for Science and Technology (FCT) under the projects Food4DIAB (EXPL/BAA-AGR/1382/2021); UIDB/50006/2020; UIDP/50006/2020 and by AgriFood XXI I&D&I project (NORTE-01-0145-FEDER-000041) co-financed by European Regional Development Fund (ERDF), through the NORTE 2020 (Programa Operacional Regional do Norte 2014/2020).

Data Availability Statement: Data is contained within the article.

Conflicts of Interest: The authors declare no conflict of interest.

### References

- Hönow, R.; Hesse, A. Comparison of Extraction Methods for the Determination of Soluble and Total Oxalate in Foods by HPLC-Enzyme-Reactor. *Food Chem.* 2002, 78, 511–521. [CrossRef]
- Li, P.; Liu, C.; Luo, Y.; Shi, H.; Li, Q.; Pinchu, C.; Li, X.; Yang, J.; Fan, W. Oxalate in Plants: Metabolism, Function, Regulation, and Application. J. Agric. Food Chem. 2022, 70, 16037–16049. [CrossRef] [PubMed]
- Petroski, W.; Minich, D.M. Is There Such a Thing as "Anti-Nutrients"? A Narrative Review of Perceived Problematic Plant Compounds. *Nutrients* 2020, 12, 2929. [CrossRef] [PubMed]
- 4. Okombo, J.; Liebman, M. Oxalate Content of Selected Breads and Crackers. J. Food Compos. Anal. 2010, 23, 118–121. [CrossRef]
- Al-Wahsh, I.A.; Wu, Y.; Liebman, M. A Comparison of Two Extraction Methods for Food Oxalate Assessment. J. Food Res. 2012, 1, 233. [CrossRef]
- Ruan, Q.-Y.; Zheng, X.-Q.; Chen, B.-L.; Xiao, Y.; Peng, X.-X.; Leung, D.W.M.; Liu, E.-E. Determination of Total Oxalate Contents of a Great Variety of Foods Commonly Available in Southern China Using an Oxalate Oxidase Prepared from Wheat Bran. J. Food Compos. Anal. 2013, 32, 6–11. [CrossRef]
- Mou, B. Evaluation of Oxalate Concentration in the U.S. Spinach Germplasm Collection. *HortScience* 2008, 43, 1690–1693. [CrossRef]
- Kasidas, G.P.; Rose, G.A. Oxalate Content of Some Commond Foods: Determination by an Enzymatic Method. *J. Hum. Nutr.* 1980, 34, 255–266. [CrossRef]
- 9. Milardović, S.; Grabarić, Z.; Rumenjak, V.; Jukić, M. Rapid Determination of Oxalate by an Amperometric Oxalate Oxidase-Based Electrode. *Electroanalysis* 2000, *12*, 1051–1058. [CrossRef]
- Bohn, T.; Davidsson, L.; Walczyk, T.; Hurrell, R.F. Fractional Magnesium Absorption Is Significantly Lower in Human Subjects from a Meal Served with an Oxalate-Rich Vegetable, Spinach, as Compared with a Meal Served with Kale, a Vegetable with a Low Oxalate Content. Br. J. Nutr. 2004, 91, 601–606. [CrossRef]
- Chai, W.; Liebman, M. Effect of Different Cooking Methods on Vegetable Oxalate Content. J. Agric. Food Chem. 2005, 53, 3027–3030. [CrossRef]
- 12. Pérez-Ruiz, T.; Martínez-Lozano, C.; Tomás, V.; Casajús, R. Flow Injection Spectrofluorimetric Determination of Oxalate Based on Its Enhancing Effect on the Oxidation of Rhodamine B by Dichromate. *Analyst* **1995**, *120*, 2111–2114. [CrossRef]

- Shaidarova, L.G.; Leksina, Y.A.; Chelnokova, I.A.; Gedmina, A.V.; Budnikov, H.C. Flow Injection Amperometric Determination of Oxalic Acid on Graphite Electrode Covered by Nafion Film with Included Palladium Nanoparticles. *Turk. Online J. Des. Art Commun.* 2017, 7, 1879–1886.
- Pundir, C.S.; Chauhan, N.; Rajneesh; Verma, M.; Ravi, A. Novel Amperometric Biosensor for Oxalate Determination Using Multi-Walled Carbon Nanotube-Gold Nanoparticle Composite. *Sens. Actuators B Chem.* 2011, 155, 796–803. [CrossRef]
- 15. Šljukić, B.; Baron, R.; Compton, R.G. Electrochemical Determination of Oxalate at Pyrolytic Graphite Electrodes. *Electroanalysis* **2007**, *19*, 918–922. [CrossRef]
- 16. Fakhari, A.R.; Rafiee, B.; Ahmar, H.; Bagheri, A. Electrocatalytic Determination of Oxalic Acid by TiO2 Nanoparticles/Multiwalled Carbon Nanotubes Modified Electrode. *Anal. Methods* **2012**, *4*, 3314–3319. [CrossRef]
- 17. Holmes, R.P.; Kennedy, M. Estimation of the Oxalate Content of Foods and Daily Oxalate Intake. *Kidney Int.* **2000**, *57*, 1662–1667. [CrossRef]
- Merusi, C.; Corradini, C.; Cavazza, A.; Borromei, C.; Salvadeo, P. Determination of Nitrates, Nitrites and Oxalates in Food Products by Capillary Electrophoresis with PH-Dependent Electroosmotic Flow Reversal. *Food Chem.* 2010, 120, 615–620. [CrossRef]
- 19. Sotomayor, M.D.P.T.; Raimundo, I.M., Jr.; Oliveira Neto, G.O.; Kubota, L.T. Bi-Enzymatic Optode Detection System for Oxalate Determination Based on a Natural Source of Enzyme. *Anal. Chim. Acta* **2001**, *447*, 33–40. [CrossRef]
- 20. Ojinnaka, M.C.; Ebinyasi, C.S.; Ihemeje, A.; Okorie, S.U. Nutritional Evaluation of Complementary Food Gruels Formulated from Blends of Soybean Flour and Ginger Modified Cocoyam Starch. *Adv. J. Food Sci. Technol.* **2013**, *5*, 1325–1330. [CrossRef]
- Adeniyi, S.A.; Orjiekwe, C.L.; Ehiagbonare, J.E. Determination of Alkaloids and Oxalates in Some Selected Food Samples in Nigeria. *Afr. J. Biotechnol.* 2009, *8*, 110–112.
- 22. Gouveia, C.S.S.; Ganança, J.F.T.; Lebot, V.; de Carvalho, M.Â.A.P. Quantitation of Oxalates in Corms and Shoots of Colocasia Esculenta (L.) Schott under Drought Conditions. *Acta Physiol. Plant* **2018**, *40*, 214. [CrossRef]
- 23. Ikese, O.; Ubwa, S.; Sunday, A.; Lenka, J.; Inalegwu, J.; Ubwa, S.; Ocheje, M.; Inegedu, A. Proximate Composition, Antinutrients and Some Functional Properties of a Potential Infant Food Made from Wheat and Groundnut. *Int. J. Food Sci. Nutr.* **2016**, *1*, 59–63.
- 24. Iwuozor, K.O. Qualitative and Quantitative Determination of Anti-Nutritional Factors of Five Wine Samples. *Adv. J. Chem.* 2019, 2, 136–146. [CrossRef]
- 25. Baker, C.J.L. The Determination of Oxalates in Fresh Plant Material. Analyst 1952, 77, 340–344. [CrossRef]
- 26. Zhang, Y.; Li, Y.; Wei, J.; Sun, M.; Tian, Y.; Li, Z. Effects of Nitrogen and Calcium Nutrition on Oxalate Contents, Forms, and Distribution in Spinach. *J. Plant Nutr.* **2009**, *32*, 2123–2139. [CrossRef]
- 27. Zhang, Y.; Lin, X.; Zhang, Y.; Shao, J.Z.; Du, S. Effects of Nitrogen Levels and Nitrate/Ammonium Ratios on Oxalate Concentrations of Different Forms in Edible Parts of Spinach. J. Plant Nutr. 2005, 28, 2011–2025. [CrossRef]
- 28. Ohkawa, H. Gas Chromatographic Determination of Oxalic Acid in Foods. J. Assoc. Off. Anal. Chem. 1985, 68, 108–111. [CrossRef]
- 29. Dong, M.W. *Modern HPLC for Practicing Scientists*; Wiley-Interscience: Hoboken, NJ, USA, 2006; ISBN 9780471727897.
- 30. Altunay, N.; Gürkan, R. A Simple, Low-Cost, and Useful Preconcentration Method for Quantification of Soluble, Insoluble, and Total Oxalate in Selected Vegetables Through Spectrophotometry. *Food Anal. Methods* **2016**, *9*, 950–965. [CrossRef]
- Ensafi, A.A.; Abbasi, S.; Rezaei, B. Kinetic Spectrophotometric Method for the Determination of Oxalic Acid by Its Catalytic Effect on the Oxidation of Safranine by Dichromate. Spectrochim. Acta Part A 2001, 57, 1833–1838. [CrossRef]
- 32. Savage, G.P.; Vanhanen, L.; Mason, S.M.; Ross, A.B. Effect of Cooking on the Soluble and Insoluble Oxalate Content of Some New Zealand Foods. *J. Food Compos. Anal.* 2000, 13, 201–206. [CrossRef]
- Radek, M.; Savage, G.P. Oxalates in Some Indian Green Leafy Vegetables. Int. J. Food Sci. Nutr. 2008, 59, 246–260. [CrossRef] [PubMed]
- Savage, G.P.; Mårtensson, L. Comparison of the Estimates of the Oxalate Content of Taro Leaves and Corms and a Selection of Indian Vegetables Following Hot Water, Hot Acid and in Vitro Extraction Methods. J. Food Compos. Anal. 2010, 23, 113–117. [CrossRef]
- 35. Wilson, C.W.; Shaw, P.E.; Knight, R.J. Analysis of Oxalic Acid in Carambola (Averrhoa Carambola L.) and Spinach by High-Performance Liquid Chromatography. *J. Agric. Food Chem.* **1982**, *30*, 1106–1108. [CrossRef]
- 36. Akhtar, M.S.; Israr, B.; Bhatty, N.; Ali, A. Effect of Cooking on Soluble and Insoluble Oxalate Contents in Selected Pakistani Vegetables and Beans. *Int. J. Food Prop.* **2011**, *14*, 241–249. [CrossRef]
- Ensafi, A.A.; Kazemzadeh, A. Flow Injection Spectrophotometric Determination of Ultra Trace Amounts of Oxalic Acid. *Fresenius J. Anal. Chem.* 2000, 367, 590–592. [CrossRef] [PubMed]
- Jiang, Z.-L.; Zhao, M.-X.; Liao, L.-X. Catalytic Spectrophotometric Methods for the Determination of Oxalic Acid. Anal. Chim. Acta 1996, 320, 139–143. [CrossRef]
- Liu, X.X.; Zhou, K.; Hu, Y.; Jin, R.; Lu, L.L.; Jin, C.W.; Lin, X.Y. Oxalate Synthesis in Leaves Is Associated with Root Uptake of Nitrate and Its Assimilation in Spinach (*Spinacia oleracea* L.) Plants. J. Sci. Food Agric. 2014, 95, 2105–2116. [CrossRef]
- Kim, D.-J.; Kim, H.; Kim, M.; Lee, J. Analysis of Oxalic Acid of Various Vegetables Consumed in Korea. Food Sci. Biotechnol. 2007, 16, 650–654.
- 41. Morosanova, M.A.; Samodelov, Z.V.; Morosanova, E.I. Determination of Food Oxalates Using Silica–Titania Xerogel Modified with Eriochrome Cyanine R. *Sensors* **2018**, *18*, 864. [CrossRef]

- 42. Fatoki, O.S. Determination of Oxalic Acid in Vegetables. In *Modern Methods of Plant Analysis: Vegetables and Vegetable Products;* Springer Berlin: Heidelberg, Germany, 1994; Volume 16, pp. 161–167.
- 43. Singh, P.P.; Kothari, L.K.; Sharma, D.C.; Saxena, S.N. Nutritional Value of Foods in Relation to Their Oxalic Acid Content. *Am. J. Clin. Nutr.* **1972**, *25*, 1147–1152. [CrossRef] [PubMed]
- Natesh, H.N.; Abbey, L.; Asiedu, S.K. An Overview of Nutritional and Anti Nutritional Factors in Green Leafy Vegetables. *Hortic.* Int. J. 2017, 1, 58–65. [CrossRef]
- 45. Huynh, N.K.; Nguyen, D.H.M.; Nguyen, H.V.H. Effects of Processing on Oxalate Contents in Plant Foods: A Review. J. Food Compos. Anal. 2022, 112, 104685. [CrossRef]
- 46. Wang, Z.; Zhang, Y.; Zhang, J.; Deng, Q.; Liang, H. Recent Advances on the Mechanisms of Kidney Stone Formation (Review). *Int. J. Mol. Med.* **2021**, *48*, 149. [CrossRef] [PubMed]
- 47. Massey, L.K. Food Oxalate: Factors Affecting Measurement, Biological Variation, and Bioavailability. J. Am. Diet. Assoc. 2007, 107, 1191–1194. [CrossRef]
- 48. Witting, C.; Langman, C.B.; Assimos, D.; Baum, M.A.; Kausz, A.; Milliner, D.; Tasian, G.; Worcester, E.; Allain, M.; West, M.; et al. Pathophysiology and Treatment of Enteric Hyperoxaluria. *Clin. J. Am. Soc. Nephrol.* **2021**, *16*, 487–495. [CrossRef]
- Sharma, S.; Rao, R.N.; Pani, K.C.; Paul, P. Bone Marrow Oxalosis: An Unusual Cause of Cytopenia in End-Stage Renal Disease; Report of Two Cases. *Indian J. Pathol. Microbiol.* 2018, 61, 268–270. [CrossRef]
- 50. Fogo, A.B.; Lusco, M.A.; Najafian, B.; Alpers, C.E. AJKD Atlas of Renal Pathology: Oxalosis. *Am. J. Kidney Dis.* **2017**, *69*, e13–e14. [CrossRef]
- 51. Noonan, S.C.; Savage, G.P. Oxalate Content of Foods and Its Effect on Humans. Asia Pac. J. Clin. Nutr. 1999, 8, 64–74. [CrossRef]
- 52. López-Moreno, M.; Garcés-Rimón, M.; Miguel, M. Antinutrients: Lectins, Goitrogens, Phytates and Oxalates, Friends or Foe? J. Funct. Foods 2022, 89, 104938. [CrossRef]
- Karamad, D.; Khosravi-Darani, K.; Hosseini, H.; Tavasoli, S. Analytical Procedures and Methods Validation for Oxalate Content Estimation. *Biointerface Res. Appl. Chem.* 2019, *9*, 4305–4310. [CrossRef] [PubMed]
- 54. Chamjangali, M.A.; Keley, V.; Bagherian, G. Kinetic Spectrophotometric Method for the Determination of Trace Amounts of Oxalate by an Activation Effect. *Anal. Sci.* **2006**, *22*, 333–336. [CrossRef] [PubMed]
- Libert, B.O. Rapid Determination of Oxalic Acid by Reversed-Phase High-Performance Liquid Chromatography. J. Chromatogr. 1981, 210, 540–543. [CrossRef]
- Arias-Rico, J.; Macías-León, F.J.; Alanís-García, E.; Cruz-Cansino, N.d.S.; Jaramillo-Morales, O.A.; Barrera-Gálvez, R.; Ramírez-Moreno, E. Study of Edible Plants: Effects of Boiling on Nutritional, Antioxidant, and Physicochemical Properties. *Foods* 2020, 9, 599. [CrossRef]
- 57. Nemzer, B.; Al-Taher, F.; Abshiru, N. Extraction and Natural Bioactive Molecules Characterization in Spinach, Kale and Purslane: A Comparative Study. *Molecules* 2021, *26*, 2515. [CrossRef]
- 58. Kassie, W.; Washe, A.P.; Etsay, H. Spectrophotometric Determination of Oxalic Acid in Dietary Sources Through Catalytic Titration with Hexavalent Chromium. *Food Sci. Qual. Manag.* **2019**, *83*, 30–38. [CrossRef]
- 59. Santamaria, P.; Elia, A.; Serio, F.; Todaro, E. A Survey of Nitrate and Oxalate Content in Fresh Vegetables. *J. Sci. Food Agric.* **1999**, 79, 1882–1888. [CrossRef]
- 60. Elia, A.; Santamaria, P.; Serio, F. Nitrogen Nutrition, Yield and Quality of Spinach. J. Sci. Food Agric. 1998, 76, 341–346. [CrossRef]
- 61. Judprasong, K.; Charoenkiatkul, S.; Sungpuag, P.; Vasanachitt, K.; Nakjamanong, Y. Total and Soluble Oxalate Contents in Thai Vegetables, Cereal Grains and Legume Seeds and Their Changes after Cooking. J. Food Compos. Anal. 2006, 19, 340–347. [CrossRef]
- Ombódi, A.; Kosuge, S.; Saigusa, M. Effects of Polyolefin-Coated Fertilizer on Nutritional Quality of Spinach Plants. J. Plant Nutr. 2000, 23, 1495–1504. [CrossRef]
- 63. Siener, R.; Seidler, A.; Hönow, R. Oxalate-Rich Foods. Food Sci. Technol. 2021, 41, 169–173. [CrossRef]
- 64. Siener, R.; Hönow, R.; Voss, S.; Seidler, A.; Hesse, A. Oxalate Content of Cereals and Cereal Products. J. Agric. Food Chem. 2006, 54, 3008–3011. [CrossRef] [PubMed]
- 65. Kusuma, D.S.; Vanhanen, L.P.; Savage, G.P. Evaluation of Extraction Parameters for Total Oxalate Determination in Spinach Using Design of Experiment Analysis. *J. Food Compos. Anal.* **2016**, *51*, 9–14. [CrossRef]
- Ghosh Das, S.; Savage, G.P. Oxalate Content of Indian Spinach Dishes Cooked in a Wok. J. Food Compos. Anal. 2013, 30, 125–129. [CrossRef]
- Vanhanen, L.; Savage, G. Comparison of Oxalate Contents and Recovery from Two Green Juices Prepared Using a Masticating Juicer or a High Speed Blender. NFS J. 2015, 1, 20–23. [CrossRef]
- Zulkhairi, A.M.; Razali, M.; Umikalsum, M.B.; Norfaizal, G.M.; Athirah, A.A.; Aisyah, M.N.S. Determination of Oxalates in Corms of Selected Taro (*Colocasia esculenta*) Varieties in Malaysia Using Ultra High-Performance Liquid Chromatography. *Asian J. Chem. Sci.* 2020, 7, 28–37. [CrossRef]
- Jaworska, G. Nitrates, Nitrites, and Oxalates in Products of Spinach and New Zealand Spinach: Effect of Technological Measures and Storage Time on the Level of Nitrates, Nitrites, and Oxalates in Frozen and Canned Products of Spinach and New Zealand Spinach. *Food Chem.* 2005, 93, 395–401. [CrossRef]
- 70. Charrier, M.J.S.; Savage, G.P.; Vanhanen, L. Oxalate Content and Calcium Binding Capacity of Tea and Herbal Teas. *Asia Pac. J. Clin. Nutr.* 2002, *11*, 298–301. [CrossRef]

- 71. Yusenko, E.; Polyntseva, E.; Lyzhova, A.; Kalyakina, O. Determination of Oxalate and Some Inorganic Anions in Green and Black Tea. *Proc. Latv. Acad. Sci. Sect. B* 2013, 67, 429–432. [CrossRef]
- 72. Siener, R.; Seidler, A.; Voss, S.; Hesse, A. Oxalate Content of Beverages. J. Food Compos. Anal. 2017, 63, 184–188. [CrossRef]
- 73. Kaminishi, A.; Kita, N. Seasonal Change of Nitrate and Oxalate Concentration in Relation to the Growth Rate of Spinach Cultivars. *HortScience* **2006**, *41*, 1589–1595. [CrossRef]
- 74. Ensafi, A.A.; Emadi, M. Spectrophotometric Reaction Rate Method for Determination of Oxalic Acid in Food Based on Its Enhancing Effect on the Oxidation of Pyrocathecol Violet by Dichromate. *Anal. Lett.* **2004**, *37*, 321–332. [CrossRef]
- 75. Tabe, M.; Fujimoto, T.; Nakahara, R.; Yamaguchi, T.; Fujita, Y. Spectrophotometric Determination of Oxalate Ion with N,N'-Diethyl-N,N'-[[4,4'-Dihydroxy-1,1'-Binaphthalene]-3,3'-Diyl]Bisbenzamide and Copper(II). Anal. Sci. 2007, 23, 601–604. [CrossRef] [PubMed]
- 76. Tavallali, H.; Deilamy-Rad, G.; Mosallanejad, N. Development of a New Colorimetric Chemosensor for Selective Determination of Urinary and Vegetable Oxalate Concentration through an Indicator Displacement Assay (IDA) in Aqueous Media. *Food Technol. Biotechnol.* 2018, 56, 329–336. [CrossRef]
- 77. Matsubara, C.; Yokoi, Y.; Tsuji, M.; Takamura, K. Flow Injection Analysis of Oxalate in Foods Using Titanium(IV)-Porphyrin Reagent. *Anal. Sci.* **1995**, *11*, 245–249. [CrossRef]
- Chamjangali, M.A.; Sharif-Razavian, L.; Yousefi, M.; Amin, A.H. Determination of Trace Amounts of Oxalate in Vegetable and Water Samples Using a New Kinetic-Catalytic Reaction System. *Spectrochim. Acta A Mol. Biomol. Spectrosc.* 2009, 73, 112–116. [CrossRef]
- 79. Xu, X.-Q.; Zhang, Z.-Q. Kinetic Spectrophotometric Determination of Oxalic Acid Based on the Catalytic Oxidation of Bromophenol Blue by Dichromate. *Mikrochim. Acta* 2000, 135, 169–172. [CrossRef]
- 80. Hassouna, M.E.M.; Elsuccary, S.A.A. Determination of Oxalate Based on Its Enhancing Effect on the Oxidation of Mn(II) by Periodate. *Talanta* 2002, *56*, 193–202. [CrossRef]
- Yan, Z.-Y.; Xing, G.-M.; Li, Z.-X. Quantitative Determination of Oxalic Acid Using Victoria Blue B Based on a Catalytic Kinetic Spectrophotometric Method. *Microchim. Acta* 2004, 144, 199–205. [CrossRef]
- 82. Safavi, A.; Banazadeh, A.R. Catalytic Determination of Traces of Oxalic Acid in Vegetables and Water Samples Using a Novel Optode. *Food Chem.* **2007**, *105*, 1106–1111. [CrossRef]
- 83. Holloway, W.D.; Argall, M.E.; Jealous, W.T.; Lee, J.A.; Bradbury, J.H. Organic Acids and Calcium Oxalate in Tropical Root Crops. J. Agric. Food Chem. 1989, 37, 337–341. [CrossRef]
- 84. Buslig, B.S.; Wilson, C.W.; Shaw, P.E. High-Performance Liquid Chromatographic Separation of Carboxylic Acids with Anion-Exchange and Reverse-Phase Columns. *J. Agric. Food Chem.* **1982**, *30*, 342–345. [CrossRef]
- 85. Glod, B.K. Ion Exclusion Chromatography: Parameters Influencing Retention. *Neurochem. Res.* **1997**, *22*, 1237–1248. [CrossRef] [PubMed]
- 86. Swartz, M. HPLC Detectors: A Brief Review. J. Liq. Chromatogr. Relat. Technol. 2010, 33, 1130–1150. [CrossRef]
- 87. Hönow, R.; Bongartz, D.; Hesse, A. An Improved HPLC-Enzyme-Reactor Method for the Determination of Oxalic Acid in Complex Matrices. *Clin. Chim. Acta* **1997**, *261*, 131–139. [CrossRef] [PubMed]
- 88. Farré, M.; Xirgu, J.; Salgado, A.; Peracaula, R.; Reig, R.; Sanz, P. Fatal Oxalic Acid Poisoning from Sorrel Soup. *Lancet.* **1989**, *334*, 1524. [CrossRef]
- 89. Bargagli, M.; Tio, M.C.; Waikar, S.S.; Ferraro, P.M. Dietary Oxalate Intake and Kidney Outcomes. *Nutrients* **2020**, *12*, 2673. [CrossRef]
- 90. Taylor, E.N.; Curhan, G.C. Oxalate Intake and the Risk for Nephrolithiasis. J. Am. Soc. Nephrol. 2007, 18, 2198–2204. [CrossRef]
- 91. Mitchell, T.; Kumar, P.; Reddy, T.; Wood, K.D.; Knight, J.; Assimos, D.G.; Holmes, R.P. Dietary Oxalate and Kidney Stone Formation. *Am. J. Physiol.-Ren. Physiol.* **2019**, *316*, 409–413. [CrossRef]
- 92. Huang, Y.; Zhang, Y.H.; Chi, Z.P.; Huang, R.; Huang, H.; Liu, G.Y.; Zhang, Y.F.; Yang, H.S.; Lin, J.H.; Yang, T.H.; et al. The Handling of Oxalate in the Body and the Origin of Oxalate in Calcium Oxalate Stones. *Urol. Int.* **2020**, *104*, 167–176. [CrossRef]
- Karamad, D.; Khosravi-Darani, K.; Khaneghah, A.M.; Miller, A.W. Probiotic Oxalate-Degrading Bacteria: New Insight of Environmental Variables and Expression of the Oxc and Frc Genes on Oxalate Degradation Activity. *Foods* 2022, 11, 2876. [CrossRef] [PubMed]
- Wang, B.; Wu, B.; Liu, J.; Yao, W.; Xia, D.; Li, L.; Chen, Z.; Ye, Z.; Yu, X. Analysis of Altered MicroRNA Expression Profiles in Proximal Renal Tubular Cells in Response to Calcium Oxalate Monohydrate Crystal Adhesion: Implications for Kidney Stone Disease. *PLoS ONE* 2014, 9, e0101306. [CrossRef] [PubMed]
- Ryall, R.L.; Fleming, D.E.; Doyle, I.R.; Evans, N.A.; Dean, C.J.; Marshall, V.R. Intracrystalline Proteins and the Hidden Ultrastructure of Calcium Oxalate Urinary Crystals: Implications for Kidney Stone Formation. *J. Struct. Biol.* 2001, 134, 5–14. [CrossRef] [PubMed]
- 96. Panis, V.; Tosios, K.I.; Gagari, E.; Griffin, T.J.; Damoulis, P.D. Severe Periodontitis in a Patient With Hyperoxaluria and Oxalosis: A Case Report and Review of the Literature. *J. Periodontol.* **2010**, *81*, 1497–1504. [CrossRef] [PubMed]
- 97. Bacchetta, J.; Farlay, D.; Abelin-Genevois, K.; Lebourg, L.; Cochat, P.; Boivin, G. Bone Impairment in Oxalosis: An Ultrastructural Bone Analysis. *Bone* **2015**, *81*, 161–167. [CrossRef]
- 98. D'Costa, M.R.; Winkler, N.S.; Milliner, D.S.; Norby, S.M.; Hickson, L.T.J.; Lieske, J.C. Oxalosis Associated With High-Dose Vitamin C Ingestion in a Peritoneal Dialysis Patient. *Am. J. Kidney Dis.* **2019**, *74*, 417–420. [CrossRef]

- Blackmon, J.A.; Jeffy, B.G.; Malone, J.C.; Knable, A.L., Jr. Oxalosis Involving the Skin: Case Report and Literature Review. Arch. Dermatol. 2011, 147, 1302–1305. [CrossRef]
- Lorenz, E.C.; Michet, C.J.; Milliner, D.S.; Lieske, J.C. Update on Oxalate Crystal Disease. Curr. Rheumatol. Rep. 2013, 15, 340.
   [CrossRef]
- 101. Bacchetta, J.; Boivin, G.; Cochat, P. Bone Impairment in Primary Hyperoxaluria: A Review. *Pediatr. Nephrol.* **2016**, *31*, 1–6. [CrossRef]
- 102. Bhasin, B. Primary and Secondary Hyperoxaluria: Understanding the Enigma. World J. Nephrol. 2015, 4, 235–244. [CrossRef]
- Soliman, N.R.; Effat, B.A.M.; Mehanna, N.S.; Tawfik, N.F.; Ibrahim, M.K. Activity of Probiotics from Food Origin for Oxalate Degradation. *Arch. Microbiol.* 2021, 203, 5017–5028. [CrossRef]
- 104. Ermer, T.; Nazzal, L.; Tio, M.C.; Waikar, S.; Aronson, P.S.; Knauf, F. Oxalate Homeostasis. *Nat. Rev. Nephrol.* 2023, 19, 123–138. [CrossRef]
- 105. Food and Agriculture Organization of the United Nations; World Health Organization. *Probiotics in Food: Health and Nutritional Properties and Guidelines for Evaluation*; Food and Agriculture Organization of the United Nations: Rome, Italy, 2006.
- 106. Gomathi, S.; Sasikumar, P.; Anbazhagan, K.; Sasikumar, S.; Kavitha, M.; Selvi, M.S.; Selvam, G.S. Screening of Indigenous Oxalate Degrading Lactic Acid Bacteria from Human Faeces and South Indian Fermented Foods: Assessment of Probiotic Potential. *Sci. World J.* 2014, 2014, 648059. [CrossRef]

**Disclaimer/Publisher's Note:** The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.